

# Gated Integrators and Boxcar Averagers

SR245 — Computer interface module with GPIB and RS-232



- **Eight analog I/O ports**
- **8-bit digital I/O port**
- **Two TTL I/O ports**
- **RS-232 and GPIB interfaces**
- **3500 point sample memory**
- **Simple command structure**

## SR245 Computer Interface

The SR245 Computer Interface module is a powerful tool for data acquisition. It provides both an analog and a digital interface between your computer and your experiment.

### Analog I/O

The eight analog I/O channels can be designated through software as all inputs, all outputs, or as a combination of inputs and outputs. All channels have 13 bits of resolution over the  $\pm 10.24$  VDC full-scale range, with 0.05 % accuracy.

### Digital I/O

Two front-panel digital I/O bits are provided for use as counters or triggers and can be set or read by the computer. Additionally, an 8-bit input and an 8-bit output port are available (on an internal connector) for your own custom digital interfaces.

### RS-232 and GPIB interfaces

Both RS-232 and GPIB interfaces are standard features of the SR245. Simple commands make programming easy from a variety of high-level languages—all that's necessary is the ability to send and receive ASCII strings. For example, sending "?5" instructs the module to measure the voltage on the 5<sup>th</sup> analog input BNC. Other commands allow you to record in the module's 3500 point buffer memory, ramp an analog output at a specified rate (for gate scanning), or read the contents of a digital counter.

### Ordering Information

SR245      Computer interface module

[www.gmp.ch](http://www.gmp.ch)

GMP SA Main office: Avenue des Baumettes 17

CH -1020 Renens

Tél. 021 633 21 21

Fax. 021 633 21 29

info@gmp.ch

GMP SA Büro Zürich: Dübendorfstrasse 11a

CH-8117 Fällanden

Tel. 044 825 34 00

Fax. 044 825 34 01

info@gmp.ch

**Analog Ports**

Configuration	Any number of the eight ports may be designated under program control as input ports, the rest default to output ports.
Inputs	1 MΩ impedance, ±10.24 VDC range Protected to 40 VDC 13-bit resolution (2.5 mV) 0.5 % accuracy Input offset <2.5 mV Maximum A/D rate is 2 kHz
Outputs	<1 Ω impedance Short circuit current limit is 20 mA 13-bit resolution (2.5 mV) 0.5 % accuracy Output offset <2.5 mV

**Digital Ports**

Type	Two front-panel I/O TTL bits, one 8-bit digital input port, one 8-bit latched digital output port
Front-panel inputs	Input impedances >100 kΩ Minimum pulse width is 200 ns Maximum count rate is 4 MHz Logic one >3 VDC, logic zero <0.7 VDC Inputs protected to ±10 VDC
Front-panel outputs	Can drive 50 Ω loads to TTL levels

**General**

Interfaces	IEEE-488 (GPIB) and RS-232 (110 baud to 19.2 kbaud)
Power	+24 V/60 mA, 24 V/60 mA, +12 V/20 mA, approx. 8 watts
Mechanical Warranty	Single-width standard NIM module One year parts and labor on defects in materials and workmanship

**Command List**

*Input/Output Commands*

I<n> n=0 to 8	Designates the first n analog ports as inputs, the remainder become outputs
?<n> n=1 to 8	Returns the value of the designated analog port
?B<n> n=1,2	Returns the value (0 or 1) of the designated digital port
?D	Returns the value of the 8-bit digital input port
?S	Returns the value of the status byte, and clears the status byte
C	Configures B2 as an input and resets the B2 counter
?C	Returns number of pulses occurring at B2 since the previous ?C
S<n>=<x>	Sets the analog port n (which must

	be designated as an output) to the value x (x = -10.237 to +10.237 V) n = 1 to 8
SB<n>=<m>	Designates digital bit n as output and sets its value to m (n = 1, 2 and m = 0, 1)
SB<n>=I	Designates the selected bit as an input (n = 1, 2)
SD=<n>	Sets the 8-bit digital output port to the value n (n = 0 to 255)
SM=<n>	Sets the GPIB SRQ mask to the value n (n = 0 to 255)

*Trigger Commands*

MS	Sets the synchronous mode. Responses to ? commands are returned after next trigger.
MA	Sets the asynchronous mode (default). Responses to ? commands are returned after command is received.
T<n>	Designates every n <sup>th</sup> pulse at B1 as a trigger (n = 1 to 32,767)
DT	Masks the trigger input so that no triggers are recognized
ET	Unmasks the trigger input
PB<n>	Outputs a 10 μs TTL pulse at digital port n (n = 1, 2)
P/<n>	Outputs a 10 μs TTL pulse at B2 each n <sup>th</sup> trigger (n = 1 to 255)

*Scan Commands*

SC<i>,<k>:<n>	Scans the list i..k of analog ports or digital port for n triggers. Total # of samples may not exceed 3711. (i..k=1 to 8, D)
ES	Ends the current scan immediately and resets the point sending counter
N	Sends the next point of stored scan
?N	Returns # of points scanned
A<n>,<i>	Adds n × 2.5 mV to the value of analog port 8 (must be positive) on every i <sup>th</sup> trigger (n,i = 1 to 255)
SS<i>,<k>:<n>	Scans the list i..k of analog ports or digital port for n triggers. Data is sent in a 2 byte binary format while scan is in progress.(i..k = 1 to 8, D)
X	Sends the data of a stored scan in 2 byte binary format

*Miscellaneous Commands*

MR	Master reset returns the SR245 to its default values.
W<n>	Introduces a delay of (n × 400 μs) before sending each character over the RS-232 interface (n = 0 to 255)
Z<i>,<k>	Changes the end-of-record characters sent by SR245 to those specified by the ASCII codes, i...k