Process Monitoring Systems

PPR100/200/300 — 100 amu, 200 amu and 300 amu systems



· 0.01, 0.1, 1.0 or 10 Torr inlet pressure

- 3 × 10⁻⁵ Torr-L/s flow rate
- · 2 second response time (0.1 Torr)
- · 100, 200 or 300 amu systems
- Field-replaceable electron multiplier and filament
- · 6 orders of magnitude dynamic range
- · RGA Windows software

SRS Process Monitoring Systems

The PPR Process Monitoring Systems are designed for inline process monitoring and diagnosis. Two paths are provided to the residual gas analyzer (RGA): a high conductivity path for monitoring base vacuum, and a pressure reducing path for monitoring the process at operating pressure. The pressure reducing path contains a micro-hole orifice which is designed to operate at one of the following pressures: 0.01, 0.1, 1 or 10 Torr. It reduces the sample pressure to the operating pressure of the RGA (about 10⁻⁶ Torr). This pressure drop is maintained by the pumping system which consists of a hybrid turbomolecular pump and a diaphragm pump. Both pumps are oil free and will not contaminate your process.

The inlet assembly that attaches to your process chamber is pictured above. The system includes an RGA, by-pass valve assembly and Tee, a controller, turbo pump, diaphragm pump, and Windows software program for data acquisition and control.

The software is used to operate the instrument in various modes, including analog scan, histogram mode, and pressure vs. time mode. Both Faraday cup and electron multiplier detectors are standard with the PPR system. The electron multiplier provides additional sensitivity and higher scan speeds. The PPR system is shipped completely assembled and calibrated, and is ready to attach to your vacuum process chamber.

For further details, see the specifications on the RGA systems.

www.gmp.ch

GMP SA Main office: Avenue des Baumettes 17

GMP SA Büro Zürich: Dübendorfstrasse 11a

CH -1020 Renens

Tél. 021 633 21 21

Fax. 021 633 21 29

info@gmp.ch

CH-8117 Fällanden T

Tel. 044 825 34 00

Fax. 044 825 34 01

info@gmp.ch

Performance

Gas flow $\sim 3 \times 10^{-5}$ mbar-L/s with pressure

reduction inlet active

Response time 2 seconds at 0.1 mbar inlet pressure

(scales linearly with pressure)

Start-up time 8 minutes nominal

Connections

2.75" CF flange, rotatable with Inlet

through holes

Turbo pump to controller 6 ft. cable (provided)

Turbo pump to 6 ft. flexible hose (provided),

1/4" ID × 7/16" OD backing pump

RS-232C Computer interface

(28,800 baud, 9-pin D-connector)

Software **RGA** Windows application

Pumps

High vacuum Hybrid turbomolecular/drag pump,

70 L/s, ultimate pressure

 2×10^{-9} mbar

Diaphragm pump with ultimate Backing

pressure less than 1 mbar. Protection class IP44

Cooling Requires forced air cooling

General

Dimensions

Power requirements 110 VAC @ 60 Hz, 220 VAC @

50 Hz (not field selectable), 300 W

(see sample configurations)

Weight 16 lbs. (turbo pump, by-pass

valve and Tee, RGA)


Vary with configuration

33 lbs. (diaphragm pump and controller)

One year parts and labor on defects Warranty

in material or workmanship. Pump seals and diaphragm warranted for

90 days.




Ordering Information

PPR100/1 100 amu process monitoring system PPR200/1 200 amu process monitoring system PPR300/1 300 amu process monitoring system

O100HJR 200 °C heater jacket

Replacement electron multiplier O100EM **O100RF** Replacement ThO₂/Ir filament **O100RI** Replacement ionizer kit

(w/ filament)



Geometry B (inches)

www.gmp.ch

CH -1020 Renens GMP SA Main office: Avenue des Baumettes 17 Tél. 021 633 21 21 Fax. 021 633 21 29 info@gmp.ch

GMP SA Büro Zürich: Dübendorfstrasse 11a CH-8117 Fällanden Tel. 044 825 34 00 Fax. 044 825 34 01 info@gmp.ch