



Mag-NetXTM

Magnetic Field Cancellation



Mag-NetX™

Magnetic Field Cancellation



Stainless steel struts are modular and interchangeable and made in a variety of lengths to accommodate virtually every commercial SEM.

New modular corner piece design with integrated electronics for easy installation and service


Column, floor or wall mounted



Mag-NetX™ Controller

Building upon our advanced control systems engineering and technology to actively sense and cancel building floor vibrations, we now offer Mag-NetX™, an innovative system providing active compensation of magnetic field fluctuations.

Designed both for point-of-use and OEM applications, Mag-NetX is ideal for scanning and transmission electron microscopes, electron beam lithography systems, ion beam instruments, and any tools that incorporate a charged beam. Combining Mag-NetX with TMC's advanced vibration isolation systems, we can provide the ultimate control of vibration and magnetic fields.




Features & Benefits


- Helmholtz coil pairs for maximum symmetry and uniformity
- Continuous field cancellation
- Continuous field monitoring
- Set and forget operation
- Several AC and DC cancellation modes available
- 100x field improvement (typical)
- Dynamic, 100 μs response
- Accurate field measurement
- Graphical User Interface with continuous system monitoring and analysis
- Optional feedforward compensation of line frequency and harmonics
- Optional feedforward capability for other inputs
- Optional custom field creation while suppressing disturbance
- Easy to assemble stainless steel cage, in-room wall-mount systems also available

How to order:

- ▶ Contact TMC. An Applications Engineer will configure a system for your unique requirements and provide a quotation.



Plot 1. • Magnitude of external sensor signal is measured against disturbance field strength of external excitation driving coil.
• Helmholtz Cage size 36 x 36 x 52 in. (91 x 91 x 132 cm)
• The best performance is at the system sensor location.



Plot 2. • Magnitude of external sensor signal is measured against disturbance field strength of external excitation driving coil.
• Helmholtz Cage size 60 x 60 x 84 in. (152 x 152 x 213 cm).
• 3-axial external sensor consists of 3 orthogonal coils around a 12 x 12 x 24 in. volume.
• Excitation coil positioned outside Helmholtz cage, external sensor coils positioned around system sensors.
• Due to cage dimensions, Z suppression is lower because Z-compensation field has lower uniformity than X and Y, but longer protected dimension (24 in. vs. 12 in. for X and Y).

GENERAL SPECIFICATIONS (may vary depending on configuration)


1. System Components:	Up to 3-axes orthogonal magnetic sensor, Mag-NetX controller, Up to 3 orthogonal pairs of coils	3. Mag-NetX Controller:	3 channels for X, Y, Z cancellation, 1U standard case
2. Performance:	Active Magnetic Field Cancellation Axes: X, Y, Z	Dimensions of controller:	17 x 9 x 1.75 in. (43 x 23 x 4.5 cm)
Sensor type, sensor noise	Fluxgate type, noise <10 pT/√Hz at 1 Hz	Operational modes:	After power-on: Automatic self test/calibration and switch to controlled mode in 1 minute, no user involvement required. Manual test / debug mode. True DC mode (compensating Earth magnetic field) Offset-DC mode (Earth magnetic field ignored) Track DC-shift due to microscope moving parts Quasi-DC mode (ignoring slow, >100 sec. fluctuations) Cancel-and-Create (simultaneously cancelling disturbance and creating custom field)
Max ambient DC field (sensor dependent)	Choice of 3 standard sensors: ± 70 μT, ± 100 μT, ± 250 μT	Types of control loops:	Analog feedback with digitally controlled gain, DC - 2 KHz Digital feedback Digital feedforward (cancels AC-line power frequency and harmonics without gain-stability limits of feedback)
Bandwidth	DC - 1 KHz typical, up to 2 KHz (depends on cage and sensor)	Front panel controls:	"OK" LED indicator: Green - OK, Yellow - Warning/Error LCD 2 x 20 symbols indicator: Show menu and status Bar-LED indicators: Show X, Y, Z real time strength of compensation field
Dynamic range - ability to cancel field (depends on cage and frequency)	± 10.0 μT typical, up to ± 100.0 μT 40 dB typ (55 dB max)	4 push buttons	For LCD menu access
Controlling volume vs. magnetic field flux density	50 m ³ at 10 μT RMS (depending on cage parameters) 1 m ³ at 50 μT RMS (standard controller is able to cancel Earth magnetic field with special cage construction) 10 m ³ at 50 μT (with external amplifier and special cage)	2 BNC sockets	For calibration testing/debugging
Field reduction ratio at sensor location (using typical console-mounting cage of X*Y*Z = 36 x 36 x 52 in. [91 x 91 x 132 cm])	40-50 dB [100x - 300x] (typical) in DC - 100 Hz 26 dB [20x] (typical) in 100 - 500 Hz (See Plot 1)	Interfaces:	USB socket (appears as COM port on PC): Graphical user interface for advanced tuning/testing, accepts ASCII commands and shows menu Auxiliary analog Inputs (rear DB-37): Can be used as feedforward or to create custom field GO - NO GO signal (relay): Binary, for usage as input for protected system
Field reduction ratio in a typical volume of electron microscope column: X*Y*Z = 12 x 12 x 24 in. [30 x 30 x 60 cm], using typical floor-standing cage of 60 x 60 x 84 in. (152 x 152 x 213 cm)	X, Y: 30 dB [32 x] (typical) in DC - 100 Hz X, Y: 20 dB [10 x] (typical) in 100 - 500 Hz Z: 15 dB [5.2 x] (typical) in DC - 100 Hz Z: 10 dB [3 x] (typical) in 100 - 500 Hz 0 dB at 1000 Hz (See Plot 2)	Supply voltage, power consumption	90 - 240 VAC 50/60 Hz, 500 VA max

Column-Mounted Helmholtz Coils




Column-mounted Helmholtz Coils are readily adapted to SEM columns but impractical for TEMs.

Floor-Mounted Helmholtz Coils



Floor-mounted Helmholtz Coils may be used for both SEMs and TEMs.

Helmholtz Coils on Leg Frame



Helmholtz Coils may be mounted on a TMC leg frame.

Wall-Mounted Helmholtz Coils



Wall-mounted coils are a practical alternative to column and floor-mounted coils for TEMs and SEMs installed near the center of a room.

www.gmp.ch