APPLICATION NOTE

TEMPERATURE DEPENDENCE OF PHOSPHORS VIA PHOTO- AND THERMOLUMINESCENCE

AN_P38 v.3; Mar. 17, Georgios Arnaoutakis

INTRODUCTION

 \mathcal{M}

Phosphors used in discharge lamps and light emitting diodes¹, displays² and scintillators^{3,4} are commonly characterised via photoluminescence spectroscopy. Their light emitting properties are connected to the structure and the interaction of the active ions with the lattice. Moreover, vibrations in the



FLS980 Photoluminescence Spectrometer

lattice are sensitive to environmental temperature changes. In this note, the temperature dependence of a phosphor widely used in plasma displays, Hg and Xe discharge lamps, is studied via photoluminescence and thermallyactivated luminescence or often called thermoluminescence spectroscopy.

METHODS & MATERIALS

Luminescence temperature measurements were performed using an FLS980 Photoluminescence Spectrometer equipped with double excitation and emission monochromators in L geometry, a photomultiplier tube detector (Hamamatsu, R928P) and a 450 W Xe lamp. Gratings blazed at 250 nm were used in the excitation and at 400 nm in the emission arm, with higher diffraction orders being filtered by the integrated long wave-pass filters in the FLS980.

 $BaMgAl_{10}O_{7}$: Eu phosphor (General Electric, BAM) was used without further purification treatment. The phosphor was diluted in water and deposited by evaporation on rectangular copper plates. A liquid nitrogen-cooled cryostat (Oxford Instruments, Optistat DN) was used to vary the sample temperature from 100 K to 480 K with accuracy of ± 2 K. Photoluminescence was measured under steady state excitation, while for the thermoluminescence glow curves, the sample was charged with 255 nm light for 30 s prior to the measurement. Then, illumination was blocked and various ramp rates were selected by adjusting the temperature step and stabilisation time in the spectrometer's software, F980.

RESULTS - DISCUSSION

The photoluminescence spectra of $BaMgAl_{10}O_{17}$: Eu between 120 K and 480 K are shown in Figure 1. The band with a peak around 450 nm is assigned to the $4f^{6}5d^{1}$ excited state to $4f^{7}$ ground state transition of Eu²⁺. Lower energy photoluminescence was not

detected, confirming the divalent state⁵ of europium. It can be seen that the intensity decreases monotonically as the temperature rises, associated with temperature quenching². A shift of the peak can also be seen as the temperature rises, with a maximum shift of 10 nm. Moreover, with increasing temperature, band broadening of the emission occurs with increased FWHM from approximately 40 nm to 60 nm due to reduced lattice vibrations. This photoluminescence behaviour was linked to the location of Eu ions in Beevers-Ross sites and to electron traps in the lattice⁶.



Figure 1: Photoluminescence spectra of BaMgAl₁₀O₁₇: Eu upon 255 nm excitation from 120 K to 480 K. The parameters of this measurement were, λ_{exc} =255 nm, $\Delta\lambda_{exc}$ =3 nm, $\Delta\lambda_{em}$ =1 nm, step=2 nm, dwell=0.2 s, t_{stab}=60 s.

To investigate further the role of traps in the $BaMgAl_{10}O_{17}$ lattice, the thermoluminescence of $BaMgAl_{10}O_{17}$: Eu was measured. The obtained normalised glow curves are shown in Figure 2. Linear heating rates of 30, 45 and 60 K/min were used for this measurement, without prior thermal cleaning during the heating ramps. A thermal cleaning step can be introduced by cooling the sample without illumination during this process. This way, it is ensured that any energy stored in deep traps from prior charging or environmental factors, is extracted.

A single peak can be seen at approximately 175 K, associated with emission from traps that were not thermally cleaned. The presence of the glow peak indicates the existence of traps below the band gap² of BaMgAl₁₀O₁₇. An apparent shift of the peak to higher temperatures can be seen for faster heating rates, in agreement with previous studies⁷. Furthermore, a single peak is only detected after charging at 255 nm. A structure of multiple traps is not observed at this excitation energy. On the other hand, a re-convoluted structure could be seen when charged with light above the band gap^{6,7}. This suggests that to resolve this structure, measurements at longer charging periods are required.

In addition, the technique can be readily applied to thermoluminescence excitation to confirm the origin of the emission from excited 5d states in luminescent and other phosphors¹.



EDINBURGH INSTRUMENTS



Figure 2: Normalised thermoluminescence glow curve of BaMgAl₁₀O₁₇: Eu at various heating rates. Experimental parameters were λ_{exc} =255 nm for 30 s, λ_{em} =450 nm, $\Delta\lambda_{exc}$ =10 nm, $\Delta\lambda_{em}$ =15 nm, step=5 K, dwell=0.5 s, tstab=10 s.

CONCLUSION

Photoluminescence and thermoluminescence measurements of $BaMgAl_{10}O_{17}$: Eu were presented in this application note. The temperature dependent photoluminescence was shown by band broadening and lattice relaxation, while information about the trap states were clearly shown in the thermoluminescence curves. This makes these techniques a very versatile tool for the study of photo- and thermoluminescent materials with an FLS980 Photoluminescence Spectrometer equipped with a liquid nitrogen cryostat.

REFERENCES

- Ueda, J., Dorenbos, P., Bos, A. J. J., Meijerink, A. & Tanabe, S. Insight into the Thermal Quenching Mechanism for Y₃Al₅O₁₂:Ce³⁺ through Thermoluminescence Excitation Spectroscopy. J. Phys. Chem. C 119, 25003–25008 (2015).
- Bizarri, G. & Moine, B. On the role of traps in the BaMgAl₁₀O₁₇:Eu²⁺ fluorescence mechanisms. J. Lumin. 115, 53–61 (2005).
- Bos, A. J. J., van Duijvenvoorde, R. M., van der Kolk, E., Drozdowski, W. & Dorenbos, P. Thermoluminescence excitation spectroscopy: A versatile technique to study persistent luminescence phosphors. J. Lumin. 131, 1465–1471 (2011).
- Ogiegło, J. M. et al. Luminescence and Luminescence Quenching in Gd₃(Ga,Al)₅O₁₂ Scintillators Doped with Ce³⁺. J. Phys. Chem. A 117, 2479–2484 (2013).
- Kim, Y.-I., Kim, K.-B., Jung, M.-J. & Hong, J.-S. Combined Rietveld refinement of BaMgAl₁₀O₁₇:Eu²⁺ using X-ray and neutron powder diffraction data. J. Lumin. 99, 91–100 (2002).
- Tanno, H., Zhang, S., Shinoda, T. & Kajiyama, H. Characteristics of photoluminescence, thermoluminescence and thermal degradation in Eu-doped BaMgAl₁₀O₁₇ and SrMgAl₁₀O₁₇. J. Lumin. 130, 82–86 (2010).
- Jüstel, T., Lade, H., Mayr, W., Meijerink, A. & Wiechert, D. U. Thermoluminescence spectroscopy of Eu²⁺ and Mn²⁺ doped BaMgAl₁₀O₁₇. J. Lumin. 101, 195–210 (2003).

For more information, contact:

Edinburgh Instruments

2 Bain Square, Kirkton Campus, Livingston, EH54 7DQ United Kingdom.

T: +44 (0)1506 425 300 F: +44 (0)1506 425 320

E: sales@edinst.com

W: www.edinst.com



+44 (0)1506 425 300 | sales@edinst.com | www.edinst.com

Words and logos marked with ® or ™ are registered trademarks or trademarks owned by EDINBURGH INSTRUMENTS Limited. Other brands and names mentioned herein may be the trademarks of their respective owners. Neither the whole nor any part of the information contained in, or the product described in, this document may be adapted or reproduced in any material form except with the prior written permission of the copyright holder.