



Lasers and Laser Systems

Picosecond

ADVANCED LASER TECHNOLOGIES



About Company

Background

EKSPLA focuses on the design and manufacturing of advanced lasers & systems and employs 30 years' experience as well as a close partnership with the scientific community. 76 out of the 100 top universities use EKSPLA lasers. The company is leading in the global market for scientific picosecond lasers.

Clients like CERN, NASA, ELI, Max Planck Institutes, Cambridge University and Massachusetts Institute of Technology have chosen Ekspla as their partner. For scientist who needs unique instrument for research, we provide parameter tailored laser systems that enable customer to perform complex experiments. In-house design and manufacturing ensures operative design, manufacturing and customization of new products.

Highly stable and reliable EKSPLA lasers combined with our own subsidiaries in the US, UK and China as well as more than 20 approved representative offices with properly trained laser engineers worldwide, ensure short response time and fast laser service as well as maintenance.

History

EKSPLA was founded about 30 years ago by a small team of engineers united around the idea of making the most advanced lasers in the world. EKSPLA was independent company with little money, but lots of creativity, and a deep technical understanding of lasers and how useful they could be for research and industry. From the start, the whole team had a deep mutual respect and believed in and supported each other. The first laser was sold at its first launch event, at an international exhibition in Germany. Soon after, the innovation was noticed by partners in Japan, and supply of the systems to leading universities there has been started. The concept of continuous improvement was admired and embraced, so it has become one of the key principles that apply to everything is done.

Femtosecond Lasers

Picosecond Lasers

Picosecond Tunable Systems



Contents

Femtosecond Lasers	5
FemtoLux 30 series Femtosecond Industrial Lasers	5
FemtoLux 3 series Microjoule Class Femtosecond Industrial Lasers	9
Picosecond Lasers	15
PL2210 series Diode Pumped Picosecond kHz Pulsed Nd:YAG Lasers	16
PL2230 series Diode Pumped High Energy Picosecond Nd:YAG Lasers	19
PL2250 series Flash-Lamp Pumped Picosecond Nd:YAG Lasers	23
Picosecond Tunable Systems	27
PGx01 series High Energy Broadly Tunable OPA	28
PGx11 series Transform Limited Broadly Tunable Picosecond OPA	32
PT277-XIR series Single Housing MIR Tunable Picosecond Laser	37
PT277 series Single Housing NIR-IR Range Tunable Picosecond Laser	39
PT403 series	<i>A</i> 1

Nanosecond Lasers	47
NL200 series	
Compact Q-switched DPSS Lasers	48

NL230 series High Energy Q-switched DPSS Nd:YAG Lasers	51
NL300 series Compact Flash-lamp Pumped Q-switched Nd:YAG Lasers	54
Harmonic Generators & Attenuators	

57

Tunable Wavelength Lasers 61

NT230 series High Energy Broadly Tunable DPSS Lasers	62
NT240 series Broadly Tunable kHz Pulsed DPSS Lasers	66
NT250 series Tunable Wavelength UV-NIR Range DPSS Lasers	70
NT270 series Tunable Wavelength NIR-IR Range DPSS Lasers	73
NT340 series High Energy Broadly Tunable Lasers	76

Other Ekspla Products 81

SFG spectrometer

For NL300 Series Lasers

Picosecond Vibrational Sum Frequency Generation Spectrometer 83



Applications Examples	88
Plasma Physics	89
Photolysis	90
SFG Spectroscopy	91
SHG spectroscopy & microscopy	92
Pump-Probe Spectroscopy	93
Solid-phase Photoluminescence Spectroscopy	94
Gas-phase Ion Luminescence Spectroscopy	95
Supercontinuum Generation	96
Time-Resolved Photoconductivity	97
Z-SCAN	98
LIDAR	99
LIBS	100
Terahertz Spectroscopy	101
OPCPA Seeding	102
Ultra-High Intensity	102
Photoacoustic Imaging	103

Ordering Information

104







FEMTOSECOND LASERS

FemtoLux 30 • FemtoLux 3 • UltraFlux 1k

FemtoLux 30



INNOVATIVE "DRY" COOLING SYSTEM

The FemtoLux 30 laser employs an innovative cooling system and sets new reliability standards among industrial femtosecond lasers. No additional bulky and heavy water chiller is needed.

The chiller requires periodic maintenance – cooling system draining and rinsing and water and particle filter replacement. Moreover, water leakage can cause damage to the laser head and other equipment. Instead of using water for transferring heat from a laser head, the FemtoLux 30 laser uses an innovative Direct Refrigerant Cooling method.

The refrigerant agent circulates from a PSU-integrated compressor and condenser, to a cooling plate via armored flexible lines.

The entire cooling circuit is permanently hermetically sealed and requires no maintenance.

SIMPLE & RELIABLE COOLING PLATE ATTACHMENT

The cooling plate is detachable from the laser head for more convenient laser installation.

The laser cooling equipment is integrated with the laser

power supply unit into a single 4U rack-mounted housing with a total weight of 15 kg.

PERFECT AND VERSATILE TOOL FOR MICROMACHINING

The FemtoLux 30 femtosecond laser has a tunable pulse duration from <350 fs to 1 ps and can operate in a broad AOM controlled range of pulse repetition rates from a single shot to 4 MHz.

The maximum pulse energy is more than 90 μ J operating with single pulses and can reach 250 μ J in burst mode, ensuring higher ablation rates and processing throughput for different materials.

The FemtoLux 30 beam parameters will meet the requirements of the most demanding materials and micromachining applications.

Innovative laser control electronics ensure simple control of the FemtoLux 30 laser by external controllers that could run on different platforms, be it Windows, Linux or others using REST API commands.

This makes easy integration and reduces the time and human resources required to integrate this laser into any laser micromachining equipment.

Femtosecond Industrial Lasers

FEATURES

- ▶ 30 W typical max output power
- > **90 µJ** max pulse energy
- > 250 µJ in a burst mode
- ▶ < 350 fs 1 ps
- Single shot to 4 MHz (AOM controlled)
- <0.5% RMS power long term stability over 100 hours
- ▶ M² < 1.2
- Beam circularity > 0.85
- Zero maintenance
- Dry cooling (no water used)
- PSU and cooling unit integrated into single 4U rack housing
- Easy and quick installation
- Compatible with galvo and Polygon scanners as well as PSO controllers
- ▶ 2 years of total warranty

APPLICATIONS

- LCD, LED, OLED drilling, cutting and repair
- Microelectronics manufacturing
- Glass, sapphire and ceramics micro processing
- Glass intra volume structuring
- Micro processing of different polymers and metals





Picosecond Lasers

Picosecond Tunable Systems

Nanosecond Lasers

Nanosecond Tunable Lasers

Other Ekspla Products

SPECIFICATIONS ¹⁾

Model	FemtoLux 30	
MAIN SPECIFICATIONS		
Wavelength	1030 nm	
Pulse Repetition Rate (PRR) ²⁾	200 kHz – 4 MHz	
Pulse repetition frequency (PRF) after frequency divider	PRF = PRR / N, N=1, 2, 3, , 65000; single shot	
Maximal average output power	> 27 W (typical 30 W)	
Maximal pulse energy	> 90 µJ	
Maximal total energy in a burst mode ³⁾	> 250 µJ	
Power long term stability (Std. dev.) 4)	< 0.5 %	
Pulse energy stability (Std. dev.) 5)	< 1 %	
Pulse duration (FWHM)	Tunable, < 350 fs ⁶⁾ – 1 ps	
Beam quality	M ² < 1.2 (typical < 1.1)	
Beam circularity, far field	> 0.85	
Beam divergence (full angle)	< 1 mrad	
Beam pointing thermal stability	< 20 µrad/°C	
Triggering mode	internal / external	
Pulse output control	frequency divider, pulse picker, burst mode, packet triggering, power attenuation	
Control interfaces	USB / RS232 / LAN	
Length of the umbilical cord	3 m, detachable	
Laser head cooling type	dry (direct refrigerant cooling through detachable cooling plate)	
PHYSICAL CHARACTERISTICS		
Laser head (W \times L \times H)	429 × 569 × 130 mm	
Power supply unit (W \times L \times H)	449 × 376 × 177 mm	
OPERATING REQUIREMENTS		
Mains requirements	100 – 240 V AC, single phase, 50/60 Hz	
Operating ambient temperature	18 – 27 °C	
Relative humidity	10-80 % (non-condensing)	

Air contamination level

- ¹⁾ Due to continuous improvement, all specifications are subject to change without notice. Parameters marked typical are not specifications. They are indications of typical performance and will vary with each unit we manufacture. All parameters are specified for a shortest pulse duration.
- When frequency divider is set to transmit every pulse. Fully controllable by integrated AOM.
- ³⁾ When number of pulses within a burst is set to 10 and PRR is set to a minimum value. Separation between pulses within a burst - ~20 ns.

ISO 9 (room air) or better

- ⁴⁾ Over 100 h after warm-up under constant environmental conditions.
- ⁵⁾ Under constant environmental conditions.
- At PRR > 500 kHz. At PRR < 500 kHz shortest pulse duration is < 400 fs.



SEAMLESS USER EXPERIENCE

- Easy integration.
 Remote control using REST API commands via USB, RS232 and LAN
- Reduced integration time.
 Demo electronics is available for laser control programming in advance
- Easy and quick installation.
 No water, fully disconnected laser head.
 Can be installed by the end-user.
- Easy troubleshooting.
 Integrated detectors and constant system status logging.
- No periodic maintenance required

6



PERFORMANCE



Fig 1. Typical dependence of output power and pulse energy of FemtoLux 30 laser on pulse repetition rate





Fig 5. Typical beam direction stability of FemtoLux 30 under harsh environmental conditions



Fig 2. Typical FemtoLux 30 laser (at 1030 nm) output pulse autocorrelation function







Fig 6. Average output power dependance of FemtoLux 30 laser on ambient temperature at 1030 nm



DIRECT REFRIGERANT COOLING SYSTEM

FEATURES

- Military-grade reliability
- Permanently hermetically sealed system >90,000 hour MTBF
- ▶ No maintenance
- ▶ High cooling efficiency
- >45% lower power consumption compared to water cooling equipment
- ▶ Compact and light





Simple and reliable cooling plate attachment



DRAWINGS



Fig 7. FemtoLux 30 laser head outline drawing

483

Fig 8. Power supply outline drawing

Picosecond Lasers

8



FEMTOSECOND LASERS

FemtoLux 30 • FemtoLux 3 • UltraFlux 1k

FemtoLux 3



FemtoLux 3 is a modern femtosecond fiber laser aimed for both R&D use and industrial integration. Tunable pulse duration in a range of 300 fs -5 ps, adjustable pulse repetition rate up to 10 MHz and adjustable pulse energy up to 3 µJ allows optimization of laser parameters for the desired application. These include marking and volume structuring of transparent materials, photopolymerization, biological imaging, nonlinear microscopy and many others. To expand the scope of applications even further this laser can be equipped with a second harmonics module.

With burst mode enabled, FemtoLux 3 can generate bursts of pulses with energy above 10 μ J with instant burst shape control which can significantly improve the efficiency of some processes.

Having a rigid, compact, passive air-cooled laser head and the possibility to control the laser from a wireless tablet, FemtoLux 3 can be integrated with different equipment, be it laser equipment for material micro-processing, microscopy or any other research equipment.



Microjoule Class Femtosecond Industrial Lasers

FEATURES

- ► 300 fs ... 5 ps tunable pulse duration
- Output power 3 W at 1030 nm or 1.2 W at 515 nm
- Up to 3 μJ/pulse and
 10 μJ/burst (at 1030 nm)
- Up to 1.2 µJ/pulse and 5 µJ/burst (at 515 nm)
- ▶ Excellent beam quality M² < 1.2
- Versatile laser control and syncronisation capabilities
- Up to 10 MHz pulse repetition rate
- Smart triggering for synchronous operation with polygon scanner and PSO
- Instant amplitude control
- Passive air cooling of the laser head
- ▶ 24/7 operation

APPLICATIONS

- Inner volume marking of transparent materials
- ▶ Marking and structuring
- Micromachining of brittle materials
- Photopolymerization
- Ophthalmologic surgery
- Biological Imaging
- Pumping of femtosecond OPO/OPA
- Microscopy





SPECIFICATIONS ¹⁾

1030 ± 2 nm 515 ± 1 nm			
1030 ± 2 nm 515 ± 1 nm			
1030 ± 2 nm 515 ± 1 nm			
515 ± 1 nm	1030 ± 2 nm		
	515 ± 1 nm		
< 300 fs (typical ~230 fs)			
300 fs = 5 ns			
> 3 W			
> 1.2 W			
≤ 0.5 %			
لµ 3 <			
> 1.2 µJ			
< 2 %			
1 – 10 MHz			
PRR = PRR, / N, N=1, 2, 3, , min 10 kH	lz		
via TTL input			
1 – 10 pulses			
> 10 µJ			
> 5 µJ			
via analog input			
0 – 100 % from remote control application or via a	analog input		
linear, vertical			
>1000:1			
< 1.2			
<1.0 mrad			
> 0.85			
< 30 µrad			
aperture			
2.0 ± 0.3 mm			
1.0 ± 0.2 mm			
100–240 V AC, single phase 47–63 Hz	,		
< 500 W			
15 – 30 °C			
10 – 80 % (non-condensing)			
ISO 9 (room air) or better			
air. passive			
464 x 363 v 129 mm			
620 × 363 × 129 mm			
449 × 436 × 140 mm (stand-alone) or 483 × 436 × 140 mm (19" rack mountable)			
5 m			
CLASS 4 laser product			
 ⁴⁾ At 1 MHz PRR_L under constant environmental conditions. ⁵⁾ When pulse picker is set to transmit every pulse. ⁶⁾ Pulse separation inside the burst is about 20 ns. ⁷⁾ Beam pointing stability is evaluated as a movement of the burst is table for burst. 	VOIBLE AND/OR INVOIBLE LASER ADATION AND DF GR SIM BOOGURET DO BRIET		
	$515 \pm 1 \text{ nm}$ $< 300 \text{ fs} (typical ~230 \text{ fs})$ $300 \text{ fs} - 5 \text{ ps}$ $> 3 W$ $> 1.2 W$ $< 0.5 %$ $=$		

²⁾ See typical power and energy curves for other pulse repetition rates at Fig 1., Fig 2. and Fig 4.

 $^{\scriptscriptstyle 3)}~$ At 1 MHz PRR, during 24 h of operation after warm-up under constant environmental conditions.



Nanosecond Lasers

Femtosecond Lasers

Picosecond Lasers

Picosecond Tunable Systems

Yb:Fiber 1030, 515 nm Max. 10 µJ, power 3 W, pulse 300 fs CLASS IV LASER PRODUCT

PERFORMANCE



Fig 1. Typical dependence of output power and pulse energy of FemtoLux 3 laser **at 1030 nm** when changing internal repetition rate of the laser



Fig 3. Typical FemtoLux 3 laser (at 1030 nm) output pulse autocorrelation function at 3 μ pulse energy. Calculated pulse duration is 230 fs



Fig 5. Average output power dependance on ambient temperature **at 1030 nm**



Fig 2. Typical dependence of output power and pulse energy of FemtoLux 3 laser **at 515 nm** on pulse repetition rate



Fig 4. Typical dependence of output power and pulse energy of FemtoLux 3 laser **at 1030 nm** when repetition rate is reduced by pulse picker. Internal repetition rate of the laser in this case is 1 MHz



Fig 6. Typical long term average output power stability of FemtoLux 3 laser **at 1030 nm** under constant environmental conditions

11





Fig 7. Typical M² measurement of FemtoLux 3 laser



Fig 8. Typical beam profiles along propagation axis of FemtoLux 3 series laser

REMOTE CONTROL 3:27 PM | Battery: 54% | WiFi signal: 60% | Mai | Englist Key state APPLICATION 5 Real Time Data Laser Control 515 0 i 1273 tion rate licker, kH 1273 User Manua None Disabled Home Page NaN 1 **%EKSPLA** Running (Ready) Stop

Fig 9. Example of FemtoLux 3 remote control application



129 112 8

> 304.8 362.3

300

DRAWINGS







Fig 12. Outline drawings of FemtoLux 3 laser head with second harmonic option







Fig 14. Outline drawings of FemtoLux 3 19" rack mountable control unit



Picosecond Lasers

The first EKSPLA picosecond laser has been sold on its first launch event in exhibition in Germany more than 28 years ago. Due to their excellent stability and high output parameters EKSPLA scientific picosecond lasers established their name as "Gold Standard" among scientific picosecond lasers.

Innovative design of new generation of picosecond modelocked lasers feature diodepumping-only technology, thus reducing maintenance costs and improving output parameters. Second, third, fourth and fifth (on some versions) harmonic options combined with various accessories, advanced electronics (for streak camera synchronization, phase-locked loop, synchronization of fs laser) and customization possibilities make these lasers well suited for many scientific applications, including optical parametric generator pumping, time-resolved spectroscopy, nonlinear spectroscopy, remote sensing, metrology... For customer convenience the laser can be operated from master device or personal computer through USB (VCP, ASCII commands), RS232 (ASCII commands), LAN (REST API) or RS232 (ASCII commands), LAN (REST API) depending on the system configuration or from remote control pad with backlit display that is easy to read even while wearing laser safety glasses.

SHORT SELECTION GUIDE

For Your convenience, table contains all available options and highest parameter values. Not all output specifications are available at the same time simultaneously. Please refer to the catalog page for exact specifications and available options.

Model	Max pulse energy at fundamental wavelength	Repetition rate, up to	Pumping	Pulse duration	Special feature	Page
PL2210	5 mJ at 1064 nm	1000 Hz	Diode pumped solid state	29±4 ps	kHz repetition rate	16
PL2230	40 mJ at 1064 nm	100 Hz	Diode pumped solid state	29 ± 4 ps	High pulse energy employing DPSS only technology	19
PL2250	100 mJ	20 Hz	Hybrid (DPSS master oscillator and flash-lamp pumped power amplifier)	29±4 ps	High pulse energy	23



PICOSECOND LASERS

PL2210 • PL2230 • PL2250

PL2210 SERIES



PL2210 series diode-pumped, air-cooled, mode-locked Nd:YAG lasers provide picosecond pulses at a kilohertz pulse repetition rate.

Short pulse duration, excellent pulse-to-pulse stability, superior beam quality makes PL2210 series diode pumped picosecond lasers well suited for many applications, including material processing, time-resolved spectroscopy, optical parametric generator pumping, and other tasks.

Flexible design

PL2210 series lasers offer a number of optional items that extend the capabilities of the laser. A pulse picker option allows control of the pulse repetition rate of the laser and operation in single-shot mode. The repetition rate and timing of pulses can be locked to an external RF source (with –PLL option) or other ultrafast laser system (with –FS option). The laser provides a triggering pulse for synchronization of the customer's equipment. A low jitter SYNC OUT pulse has a lead up to 500 ns that can be adjusted in ~0.25 ns steps from a PC. Up to 400 µs lead of triggering pulse is available as a PRETRIG feature that is designed to provide precise, very low jitter trigger pulses for a streak camera.

Built-in harmonic generators

Motorised switching of wavelength for PL2210A. Non-linear crystals mounted in temperature stabilized heaters are used for second, third and fourth high spectral purity harmonic generation. Diode Pumped Picosecond kHz Pulsed Nd:YAG Lasers

FEATURES

- ▶ High pulse energy at kHz rates
- Diode pumped solid state design
- Air cooled external water supply is not required (for PL2210A-1k only)
- ▶ Turn-key operation
- Low maintenance costs
- Optional streak camera triggering pulse with <10 ps rms jitter
- Remote control pad
- PC control
- Optional temperature stabilized second, third and fourth harmonic generators

APPLICATIONS

- Time resolved fluorescence (including streak camera measurements), pump-probe spectroscopy
- OPG/OPA/OPO pumping
- Remote Laser Sensing
- Other spectroscopic and nonlinear optics applications

Simple and convenient laser control

For customer convenience the laser can be operated from master device or personal computer through USB (VCP, ASCII commands), RS232 (ASCII commands), LAN (REST API) or RS232 (ASCII commands), LAN (REST API) depending on the system configuration or from remote control pad with backlit display that is easy to read even while wearing laser safety glasses.

Available models ¹⁾ Model Features

PL2210A-1k	Up to 900 μ J, 28 ps pulses at an up to 1 kHz repetition rate
PL2211A	Up to 5 mJ energy at a 1 kHz repetition rate at 28 ps pulses

¹⁾ Custom-built models with higher pulse energy are available on request.

Picosecond

Lasers

Other Ekspla Products



PL2210 SERIES

SPECIFICATIONS ¹⁾

Model	PI 2210A	DI 2211A		
at 1064 nm	0.9 ml	5 ml		
at 532 nm $^{2)}$	0.45 ml	2.5 ml		
at 355 nm ³⁾	0.35 ml	16 ml		
at 266 nm $^{4)}$	0.16 ml	1 ml		
Pulse energy stability (StdDev) 5)	U.I MJ			
at 1064 nm		0.5 %		
at 532 nm		0.8 %		
at 355 nm		1%		
at 266 nm		2 %		
Pulse duration (FWHM) ⁶⁾		29 + 4 ps		
Pulse repetition rate		1 kHz		
Triagering mode	int	ternal/external		
Typical TRIG1 OUT pulse delay ⁸⁾		-500 50 ns		
TRIGI OUT pulse iitter		< 0.1 ns rms		
Spatial mode ⁹⁾	Cla	ose to Gaussian		
Beam divergence ¹⁰⁾		<1 mrad		
Beam diameter ¹¹⁾	1.7 ± 0.3 mm	~3 mm		
Beam pointing stability ¹²⁾		< 30 μrad		
Pre-pulse contrast		> 200 : 1		
Polarization	li	near, >100 : 1		
PHYSICAL CHARACTERISTICS				
Laser head size (W \times L \times H) ¹³⁾	456 >	× 1031 × 249 mm		
Power supply size (W \times L \times H)	365 × 392 × 290 mm	550 × 600 × 550 ±3 mm (19" standard, MR-9)		
OPERATING REQUIREMENTS				
Water service	not re	quired, air cooled		
Relative humidity	20-80 9	% (non condensing)		
Ambient temperature		22 ± 2 °C		
Power requirements	100–240 V A	C, single phase 50/60 Hz		
Power consumption ¹⁴⁾	<1 kW	<1.5 kW		
Due to continuous improvement, all specifications subject to change without notice. Parameters ma	s are ⁷⁾ With respect to optical pulse. <1 rked provided optionally with PRETRIC	0 ps rms jitter is G feature.		
typical are not specifications. They are indications typical performance and will vary with each unit t manufacture. Unless stated otherwise, all specifica are measured at 1064 nm and for basic system with	or ve TRIG1 OUT lead or delay can be 0.25 ns steps in specified range. thout 9 Near field Gaussian fit is >90%.	adjusted with		
 options. For PL2210 series laser with –SH, -SH/TH, -SH/FH SH/TH/EH option. Outputs are not simultaneous. 	-SH/FH or values measured at the 1/e ² level at 1064 nm.			
 For PL2210 series laser with –TH, -SH/TH or -SH/TH/FH option. Outputs are not simultaneous. 	e ² point. ¹²⁾ Beam pointing stability is evaluat	ted from		
For PL2210 series laser with -SH/FH or -SH/TH/EH option. Outputs are not simultaneous	fluctuations of beam centroid po	fluctuations of beam centroid position in the far field.		

- ⁵⁾ Averaged from pulses, emitted during 30 sec time interval.
- ⁶⁾ Optional 80 or 22 ps ± 10% duration. Pulse energy specifications may differ from indicated here.
- ¹³⁾ 456×1233×249 mm (W×L×H) laser head size might be required for some optional configurations.
- ¹⁴⁾ At 1 kHz pulse repetition rate.

Femtosecond Lasers

Picosecond Lasers

Picosecond Tunable Systems

Nanosecond Lasers



PL2210 SERIES

OPTIONS

- PRETRIG provides low jitter pulse for streak camera triggering with lead/delay in -400...600 μs range and <10 ps rms jitter.</p>
- Option P80 provides 80 ps ± 10 % output pulse duration. Inquire for pulse energy specifications.
- Option P20 provides 22 ps ± 10 % output pulse duration. Inquire for pulse energy specifications.
- ▶ Option PC allows reduction of the pulse repetition rate of the PL2210 series laser by integer numbers. Single shot mode is also possible. In addition, the –PC option reduces the low-intensity quasi-CW background that is present at laser output at 1064 nm wavelength. Please note that the output of fundamental wavelength and harmonic will be reduced by approx. 20% with installation of the –PC option.

BEAM PROFILE

OUTLINE DRAWINGS



Fig 1. Typical PL2210 series laser near field beam profile at 1064 nm except PL2211A



Fig 2. Typical PL2211A laser near field beam profile at 1064 nm



Fig 3. Dimensions of PL2210 series laser head

ORDERING INFORMATION

Note: Laser must be connected to the mains electricity all the time. If there will be no mains electricity for longer that 1 hour then laser (system) needs warm up for a few hours before switching on.

PL2210A-SH/TH/FH-P20



gy Harmonic generator options: but SH → second harmonic

- TH \rightarrow third harmonic
- $FH \rightarrow fourth harmonic$

Other pptions: P80 \rightarrow 80

P20

PC

PLL

- \rightarrow 80 ps pulse duration option
 - \rightarrow 20 ps pulse duration option
- → pulse picker option
- → pulse repetition rate locking option



PICOSECOND LASERS

PL2210 • PL2230 • PL2250

PL2230 SERIES



Innovative design

The heart of the system is a diode pumped solid state (DPSS) master oscillator placed in a sealed monolithic block, producing high repetition rate pulse trains (87 MHz) with a low single pulse energy of several nJ. Diode pumped amplifiers are used for amplification of the pulse to 30 mJ or up to 40 mJ output. The high-gain regenerative amplifier has an amplification factor in the proximity of 10⁶. After the regenerative amplifier, the pulse is directed to a multipass power amplifier that is optimized for efficient stored energy extraction from the Nd:YAG rod, while maintaining a near Gaussian beam profile and low wavefront distortion. The output pulse energy can be adjusted in approximately 1% steps, while pulse-to-pulse energy stability remains at less than 0.5% rms at 1064 nm.

Angle-tuned KD*P and KDP crystals mounted in thermostabilised ovens are used for second, third, and fourth harmonic generation. Harmonic separators ensure the high spectral purity of each harmonic guided to different output ports. Built-in energy monitors continuously monitor output pulse energy. Data from the energy monitor can be seen on the remote keypad or on a PC monitor. The laser provides triggering pulses for the synchronisation of your equipment. The lead of the triggering pulse can be up to 500 ns and is user adjustable in ~0.25 ns steps from a personal computer. Up to 1000 µs lead of triggering pulse is available as a pretrigger feature. Precise pulse energy control, excellent short-term and long-term stability, and a 50 Hz repetition rate makes PL2230 series lasers an excellent choice for many demanding scientific applications.

Simple and convenient laser control

For customer convenience the laser can be operated from master device or personal computer through USB (VCP, ASCII commands), RS232 (ASCII commands), LAN (REST API) or RS232 (ASCII commands), LAN (REST API) depending on the system configuration or from remote control pad with backlit display that is easy to read even while wearing laser safety glasses. Diode Pumped High Energy Picosecond Nd:YAG Lasers

FEATURES

- Diode pumped power amplifier producing up to 40 mJ per pulse at 1064 nm
- Beam profile improvement using advanced beam shaping system
- Hermetically sealed DPSS master oscillator
- Diode pumped regenerative amplifier
- Air-cooled
- <30 ps pulse duration</p>
- Excellent pulse duration stability
- ▶ Up to **100 Hz** repetition rate
- Streak camera triggering pulse with <10 ps jitter
- Excellent beam pointing stability
- Thermo stabilized second, third or fourth harmonic generator options
- PC control
- Remote control via keypad

APPLICATIONS

- Time resolved fluorescence (including streak camera measurements)
- ► SFG/SHG spectroscopy
- Nonlinear spectroscopy
- Laser-induced breakdown spectroscopy
- OPG pumping
- Remote laser sensing
- Satellite ranging
- Other spectroscopic and nonlinear optics applications



PL2230 SERIES

SPECIFICATIONS ¹⁾

Model	PL2230-100	PL2231-100	PL2231-50	PL2231A-50
Pulse energy ²⁾				
at 1064 nm	3 mJ 12 mJ		30 mJ	40 mJ
at 532 nm ³⁾	1.3 mJ	5 mJ	13 mJ	18 mJ
at 355 nm 4)	0.9 mJ	3.5 mJ	9 mJ	13 mJ
at 266 nm ⁵⁾	0.3 mJ	1.2 mJ	3 mJ	5 mJ
at 213 nm ⁶		inquire		
Pulse energy stability (StdDey) 7)		• ···••]•···		
at 1064 nm	< 0.2 %		< 0.5 %	
at 532 nm	< 0.4 %		< 0.8 %	
at 355 nm	< 0.5 %		< 11 %	
at 266 nm	< 0.5 %		< 1.7 %	
at 200 mm	< 1.5 %		< 1.2 %	
Dulco duration (EN/LINA) 8)	× 1.J /0	20 + 4 m	< 1.5 /0	
Pulse duration (FWHW) *		29 ± 4 p	5	
Pulse duration stability 37		±1%		
Power drift	0	± 2 %		
Pulse repetition rate	0 – 100 Hz	100 Hz	50 Hz	50 Hz
Polarization		vertical, >99 % at	t 1064 nm	
Pre-pulse contrast	> 200) : 1 (peak-to-peak with res	spect to residual pulse	es)
Beam profile ¹¹⁾		close to Gaussian in ne	ar and far fields	
Beam divergence ¹²⁾	< 1.5 mrad		< 0.7 mrad	
Beam propagation ratio M ²	< 1.3		< 2.5	
Beam pointing stability ¹³⁾	≤ 10 µrad StdDev		≤ 20 µrad StdDev	
Typical beam diameter ¹⁴⁾	~ 2 mm	~ 6 m	ım	~ 7 mm
Optical pulse jitter				
Internal triggering regime ¹⁵⁾	<5	0 ps (StdDev) with respec	t to TRIG1 OUT pulse	
External triggering regime ¹⁶⁾		~3 ns (StdDev) with respec	t to SYNC IN pulse	
TRIG1 OUT pulse delay ¹⁷⁾		-500 50	ns	
Typical warm-up time	5 min		15 min	
		4561021240		
Laser head size (W × L × H)		456×1031×249	± 3 mm	
Electrical cabinet size (W \times L \times H)	85×170×41 ± 3 mm		471×391×147 ± 3 mm	
Umbilical length		2.5 m		
OPERATING REOUIREMENTS				
Cooling ¹⁸⁾		stand-alone of	chiller	
Room temperature		22+2 °C	•	
Relative humidity		20 – 80 % (non-co	ndensina)	
Power requirements	110-240 V AC 50/60 Hz	Single pha	so $110 - 2/0 \vee \Delta C = 5 \Delta$	50/60 Hz
Power consumption	< 0.15 kVA	Single prid	< 10 k/A	, 50/00 112
rower consumption	\$ 0.15 KVA		< 1.0 KVA	
Due to continuous improvement, all specifications are subject to change with notice. Parameters marked typical are no specifications. They are indications of typ performance and will vary with each uni	*) FWHM. Inquire for out 20 – 90 ps range. pt may differ from in pical 9 Measured over 1 i twe tomporature varia	r optional pulse durations in Pulse energy specifications dicated here. hour period when ambient tion is less than +1 %		DANGER VISIBLE AND/OR INVISIBLE LASER RADIATION AVOID ETE OF SINI EXPOSURE TO DIRECT
manufacture. Unless stated otherwise, al specifications are measured at 1064 nm a basic system without options.	and for ¹⁰⁾ Measured over 8 warm-up when ar	hours period after 20 min mbient temperature variation		Ndr4xG 106 s Sch 12kB v Abbation Ndr4xG 1064 nm, 532 nm, 355 nm, 266 nm Max 120 mJ, puge 28 ps CLASS IV LASER PRODUCT
²⁾ Outputs are not simultaneous.	11) Near fold Caucoia	 en fit is >80%		
For PL2230 series laser with –SH, -SH/TH FH or -SH/TH/FH option or –SH/TH/FH/ module.	H, -SH/ FiH at 100 Gaussie FiH fiH fiH fiH fiH fiH fiH fiH fiH fiH f	d Y-plane full angle measured at the 1/e ² level	¹⁶⁾ With respect to S ¹⁷⁾ TRIG1 OUT lead c	YNC IN pulse. or delay can be adjusted with
 4) FOR PL2230 series laser with –TH, -SH/TH -SH/TH/FH option or –SH/TH/FH/FiH mo 5) For PL2230 series laser with _SH/FH = - 	H or dule. 13) Beam pointing sta fluctuations of bea	ability is evaluated from am centroid position in the	U.25 ns steps in s ¹⁸⁾ Air cooled. Adequ should be provide	pecified range. Jate room air conditioning ed
TH/FH option or -SH/FH/FH/FH/FH module	far field. ¹⁴⁾ Beam diameter is	measured at 1064 nm at the		
 FOR PL2230 series laser with –SH/TH/FH/ module. Averaged from pulses, amitted during 20 	¹⁵⁾ With respect to TI	RIG1 OUT pulse. <10 ps jitter	If laser is optimi parametrical ge energy may be o	sed for pumping nerator, maximum output different than specified for
 Averaged from pulses, emitted during 30 time interval 	is provided option	nally with PRETRIG feature.	stand alone app	lication.

⁷⁾ Averaged from pulses, emitted during 30 sec time interval.

20



Nanosecond Tunable Lasers

PL2230 SERIES

Custom products, tailored for specific applications ¹⁾

Model	PL2231C-20 (inquire)		
Pulse energy ²⁾			
at 1064 nm	140 mJ		
at 532 nm ³⁾	60 mJ		
at 355 nm ⁴⁾	35 mJ		
at 266 nm ⁵⁾	15 mJ		
Pulse duration (FWHM) ⁶⁾	80 ps ± 10 %		
Pulse repetition rate	20 Hz		

¹⁾ Due to continuous improvement, all specifications are subject to change without notice. Parameters marked typical are not specifications. They are indications of typical performance and will vary with each unit we manufacture. Unless stated otherwise, all specifications are measured at 1064 nm and for basic system without options. Specifications for model PL2231C are preliminary and should be confirmed against quotation and purchase order.

- For PL2230 series laser with –SH, -SH/TH, -SH/FH or -SH/TH/FH option or –SH/TH/FH/FiH module.
 For PL2230 series laser with –TH, -SH/TH or -SH/TH/FH option or
- -SH/TH/FH/FH module.
- $^{\rm 5)}\,$ For PL2230 series laser with -SH/FH or -SH/TH/FH option or –SH/TH/FH/FiH module.
- $^{\rm 6)}\,$ FWHM. Inquire for optional pulse durations in 20 90 ps range. Pulse energy specifications may differ from indicated here.

OPTIONS

2) Outputs are not simultaneous.

Option P20 provides 20 ps ±10% output pulse duration. Pulse energies are ~ 30 % lower in comparison to the 28 ps pulse duration version. See table below for pulse energy specifications:

Model	PL2231-50	PL2231A-50
1064 nm	23 mJ	28 mJ
532 nm	9 mJ	13 mJ
355 nm	6 mJ	9 mJ
266 nm	2 mJ	4 mJ

- Option P80 provides 80 ps ± 10% output pulse duration. Pulse energy specifications are same as those of 28 ps lasers.
- Option PLL allows locking the master oscillator pulse train repetition rate to an external RF generator, enabling precise external triggering with low jitter. Inquire for more information.

BEAM PROFILE



Fig 1. Typical near field output beam profile of PL2230 model laser

21



PICOSECOND LASERS

PL2230 SERIES

OUTLINE DRAWINGS





ORDERING INFORMATION

Note: Laser must be connected to the mains electricity all the time. If there will be no mains electricity for longer that 1 hour then laser (system) needs warm up for a few hours before switching on.

PL2231-50-SH/TH/FH-P20

M	odel
Pulse rep	petition

rate in Hz

- Harmonic generator options: SH \rightarrow second harmonic TH \rightarrow third harmonic
- FH \rightarrow fourth harmonic

Other options:

- P20 \rightarrow 20±2 ps pulse duration option
- $\begin{array}{rrr} \mathsf{P80} & \rightarrow & 80 \text{ ps pulse duration} \\ & & \text{option} \end{array}$
- PLL → pulse repetition rate locking option



PL2210 • PL2230 • PL2250

PL2250 SERIES



PL2250 series lasers cost-effective design improves laser reliability and reduces running and maintenance costs.

Innovative design

The heart of the system is a diode pumped solid state (DPSS) master oscillator placed in a hermetically sealed monolithic block. The flashlamp pumped regenerative amplifier is replaced by an innovative diode pumped regenerative amplifier. Diode pumping results in negligible thermal lensing, which allows operation of the regenerative amplifier at variable repetition rates, as well as improved long-term stability and maintenance-free operation.

The optimized multiple-pass power amplifier is flashlamp pumped and is optimized for efficient amplification of pulse while maintaining a near Gaussian beam profile and low wavefront distortion. The output pulse energy can be adjusted in approximately 1% steps, at the same time as pulse-to-pulse energy stability remains less than 0.8% rms at 1064 nm.

Angle-tuned KD*P and KDP crystals mounted in thermostabilised ovens are used for second, third and fourth harmonic generation. Harmonic separators ensure the high spectral purity of each harmonic directed to different output ports.

Built-in energy monitors continuously monitor output pulse energy. Data from the energy monitor can be seen on the remote keypad or PC monitor. The laser provides several triggering pulses for synchronization of the customer's equipment. The lead or delay of the triggering pulse can be adjusted in 0.25 ns steps from the control pad or PC. Up to 1000 µs lead of triggering pulse is available as a pretrigger feature.

Precise pulse energy control, excellent short-term and long-term stability, and up to 20 Hz repetition rate makes PL2250 series lasers an excellent choice for many demanding scientific applications.

Simple and convenient laser control

For customer convenience the laser can be operated from master device or personal computer through USB (VCP, ASCII commands), RS232 (ASCII commands), LAN (REST API) or RS232 (ASCII commands), LAN (REST API) depending on the system configuration or from remote control pad with backlit display that is easy to read even while wearing laser safety glasses.

Flash-Lamp Pumped Picosecond Nd:YAG Lasers

FEATURES

- Hermetically sealed DPSS master oscillator
- Diode pumped regenerative amplifier
- Flashlamp pumped power amplifier producing up to 100 mJ per pulse at 1064 nm
- 30 ps pulse duration (20 ps optional)
- Excellent pulse duration stability
- Up to 20 Hz repetition rate
- Streak camera triggering pulse with <10 ps jitter
- Excellent beam pointing stability
- Thermo-stabilized second, third, fourth and fifth harmonic generator options
- ▶ PC control
- ▶ Remote control via keypad

APPLICATIONS

- Time resolved fluorescence (including streak camera measurements)
- ► SFG/SHG spectroscopy
- Nonlinear spectroscopy
- Laser-induced breakdown spectroscopy
- OPG pumping
- Remote laser sensing
- ► Satellite ranging
- Other spectroscopic and nonlinear optics experiments

Nanosecond Tunable Lasers

^picosecond Lasers

Picosecond Tunable Systems

Nanosecond Lasers



PICOSECOND LASERS

PL2250 SERIES

SPECIFICATIONS ¹⁾

Model	PL2251A	PL2251B	3 PL2251C			
Pulse energy						
at 1064 nm	50 mJ ²⁾	80 mJ ²⁾		100 mJ		
at 532 nm ³⁾	25 mJ	40 mJ		50 mJ		
at 355 nm 4)	15 mJ	24 mJ		30 mJ		
at 266 nm ⁵⁾	7 mJ	10 mJ		12 mJ		
at 213 nm ⁶⁾		inquire				
Pulse energy stability, (StdDev.) 7)						
at 1064 nm		< 0.8 %				
at 532 nm	<1.0 %					
at 355 nm		< 1.1 %				
at 266 nm		< 1.2 %				
Pulse duration (FWHM) ⁸⁾		29 ± 4 ps				
Pulse duration stability ⁹⁾		± 1.0 ps				
Repetition rate	20 o	or 10 Hz		10 Hz		
Polarization		linear, vertical, >9	99 %			
Pre-pulse contrast	>200 [.] 1 (r	peak-to-peak with respec	ct to residu	al pulses)		
Ontical pulse iitter	20011 (internal / extern	nal			
Internal triggering regime ¹⁰	< 50 ps	s (StdDev) with respect to	n TRIG1 OU	Tinulse		
External triggering regime ¹¹	~3 n	s (StdDev) with respect t		nulse		
SYNC OUT pulse litter ¹⁰	511	-500 50 ns	o orrite int	pulse		
SYNC OUT pulse delay ¹²⁾	-500 50 ns					
Beam divergence ¹³⁾	< 0.5 mrad					
Beam pointing stability ¹⁴⁾						
Beam diameter ¹⁵	~ 8 mm	~10 mm		~12 mm		
Typical warm-up time						
PHYSICAL CHARACTERISTICS						
Laser head size (W \times L \times H)	456×1233×249 mn 456×1031×249 m	n ±3 mm (for PL2251A, B 1m ±3 mm (for PL2251A,	B with harm B models v	onic and C models) vithout harmonic)		
Electric cabinet size (W \times L \times H)	550	×600×550 ±3 mm (19" s	tandard, M	R-9)		
Umbilical length		2.5 m				
OPERATING REQUIREMENTS						
Water consumption (may 20° C)	water cooled	water consumption (ma	ax 20 °C) <	8 l/min 2 har		
Boom temperature	water cooled	22 + 2 °C	IX. 20 C), <			
Relative humidity		20_80 % (pop-cond	ensing)			
Power requirements ¹⁶	sino	ule phase 200–240 V AC	16 A 50/6	0 H 7		
Power ¹⁷	~ 15 k\/A	2 5 k/A	, 10 A, 30/0	< 25 W/A		
Fower	< 1.3 KVA	< 2.3 KVA		< 2.3 KVA		
¹⁾ Due to continuous improvement, all specifications are subject to change without notice. Parameters marked typical are not specifications. They are indications of typical performance and will vary with each unit we manufacture. Unless stated otherwise, all specifications are measured at 1064 nm and for basic system without options.	 Averaged from pulses, emitted during 30 sec time interval. FWHM. Inquire for optional pulse durations in 20 – 90 ps range. Pulse energy specifications may differ from indicated here. Measured over 1 hour period when ambient temperature variation is less than ±1 °C 			VISIE ANO/OR INVISIE LASER RADATI AVOID PT CA SAN DPOSUME TO DARCT RELETION OF SANTERD ADARTON NAVING TOKA SZEM TS SANT, SZE AN STAN NAVING TOKA SZE MI SZANT, SZE AN SZANT MAK 100 mi, pulse 30 go CLASY U VLASH RODUCT		
 PL2251B-20 has 70 mJ at 1064 nm output energy. Inquire for these energies at other wavelengths. 	gy. ¹⁰ With respect to TRIG1 OUT pulse. <10 ps jitter s. is provided optionally with PRETRIG feature.					
For -SH option. Outputs are not simultaneous. Please inquire for pulse energies at other wavelengths.	 With respect to SYNC TRIG1 OUT lead or de 0.25 ns steps in specif 	IN pulse. lay can be adjusted with fied range.	¹⁶⁾ Three p require	ohase 208 or 380 VAC mains are ed for 50 Hz versions.		
 For -TH option. Outputs are not simultaneous. Please inquire for pulse energies at other wavelengths. 	 ¹³⁾ Average of X- and Y-r divergence values me at 1064 nm. 	blane full angle asured at the 1/e² level	¹⁷⁾ For 10	Hz version.		
⁵⁾ For -FH option. Outputs are not simultaneous. Please inquire for pulse energies at other wavelengths.	¹⁴⁾ Beam pointing stabilit fluctuations of beam of far field.	y is evaluated from centroid position in the	lf laser param	r is optimised for pumping etrical generator, maximum t energy may be different than		
⁶⁾ For PL2250 series laser with custom -FiH option.	ion. ¹⁵⁾ Beam diameter is measured at 1064 nm at the 1/e ² point			output energy may be different than specified for stand alone application.		

Other Ekspla Products

*****EKSPLA

PL2250 SERIES

OPTIONS

▶ Option P20 provides 20 ps ± 10% output pulse duration. Pulse energies are 30% lower in comparison to the 30 ps pulse duration version. Linewidth <2 cm⁻¹ at 1064 nm. See table below for pulse energy specifications:

Model	PL2251A-10	PL2251B-10	PL2251C -10
1064 nm	35 mJ	60 mJ	80 mJ
532 nm	17 mJ	30 mJ	40 mJ
355 nm	12 mJ	18 mJ	24 mJ
266 nm	5 mJ	8 mJ	10 mJ

▶ Option P80 provides 80 ps ±10% output pulse duration. Pulse energy specifications as below:

Model	PL2251A	PL2251B	PL2251C
Pulse energy at 1064 nm	70 mJ	100 mJ	160 mJ

BEAM PROFILE



Fig 1. Typical near field output beam profile of PL2250 series laser

OUTLINE DRAWINGS



ORDERING INFORMATION

Note: Laser must be connected to the mains electricity all the time. If there will be no mains electricity for longer that 1 hour then laser (system) needs warm up for a few hours before switching on.

PL2251A-50-SH/TH/FH-P20



Pulse repetition rate in Hz

Harmonic generator options: $SH \rightarrow$ second harmonic → third harmonic ΤH $FH \rightarrow fourth harmonic$

Other options:

- P20 20 ps pulse duration option \rightarrow P80 80 ps pulse duration option \rightarrow \rightarrow AW water-air heat exchanger
- option FS seeding option



Photo: Single Housing MIR Tunable Picosecond Laser PT277-XIR integrate a picosecond optical parametric oscillator and DPSS pump laser into a single compact housing.

Picosecond Tunable Systems

For researchers demanding wide tuning range, high conversion efficiency and narrow line-width, EKSPLA PG series optical parametric generators is an excellent choice. All models feature hands-free wavelength tuning, valuable optical components protection system as well as wide range of accessories and extension units.

Long-term experience and close cooperation with scientific institutions made it possible to create range of models, offering probably the widest tuning range: from 193 nm to 16000 nm. Versions, offering near transform limited line-width as well as operating at kHz repetition rates are available.

For customer convenience the wavelength can be set from personal computer through USB (VCP, ASCII commands), RS232 (ASCII commands), LAN (REST API) or RS232 (ASCII commands), LAN (REST API) depending on the system configuration or from remote control pad with backlit display that is easy to read even while wearing laser safety glasses. EKSPLA PL series picosecond mode-locked lasers are recommended for pumping of PG series Optical Parametric Generators. Combining together, researchers get complete tunable wavelength system, capable to assist researchers in wide range of spectroscopy applications: time-resolved pump-probe, nonlinear, infrared spectroscopy, laser-induced fluorescence.

SHORT SELECTION GUIDE

For Your convenience, table contains all available options and highest parameter values. Not all output specifications are available at the same time simultaneously. Please refer to the catalog page for exact specifications and available options.

Model	Output wavelength range	Max pulse repetition rate	Linewidth	Special feature	Page
PGx01	193 – 16 000 nm	50 Hz	< 6 cm ⁻¹	High peak power (>50 MW), ideal for non-linear spectroscopy	28
PGx11	193 – 16 000 nm	50 Hz or 1000 Hz	< 2 cm ⁻¹	Narrow linewidth (<0.8 cm ⁻¹ on some versions)	32
PT277-XIR	1400 – 16000 nm	87 MHz	< 5 cm ⁻¹	Picosecond MHz rate MIR range laser system	37
РТ277	1400 – 2050, 2200–4450 nm	87 MHz	< 2.5 cm ⁻¹	Optional intensity modulation up to 2 MHz	39
РТ403	210 – 2300 nm	1000 Hz	< 9 cm⁻¹	Pump laser and OPG integrated in 2-in-1 combo housing	41



PICOSECOND TUNABLE SYSTEMS

PGx01 • PGx11 • PT277-XIR • PT277 • PT403

PGx01 SERIES



Travelling Wave Optical Parametric Generators (TWOPG) are an excellent choice for researchers who need an ultra-fast tunable coherent light source from UV to mid IR.

Design

The units can be divided into several functional modules:

- optical parametric generator (OPG);
- diffraction grating based linewidth narrowing system (LNS);
- optical parametric amplifier (OPA);
- electronic control unit.

The purpose of the OPG module is to generate parametric superfluorescence (PS). Spectral properties of the PS are determined by the properties of a nonlinear crystal and usually vary with the generated wavelength. In order to produce narrowband radiation, the output from OPG is narrowed by LNS down to 6 cm⁻¹ and then used to seed OPA.

Output wavelength tuning is achieved by changing the angle of the nonlinear crystal(s) and grating. To ensure exceptional wavelength reproducibility, computerized control unit driven precise stepper motors rotate the nonlinear crystals and diffraction grating. Nonlinear crystal temperature stabilization ensures long-term stability of the output radiation wavelength.

In order to protect nonlinear crystals from damage, the pump pulse energy is monitored by built-in photodetectors, and the control unit produces an alert signal when pump pulse energy exceeds the preset value.

For customer convenience the laser can be operated from master device or personal computer through USB (VCP, ASCII commands), RS232 (ASCII commands), LAN (REST API) or RS232 (ASCII commands), LAN (REST API) depending on the system configuration or from remote control pad with backlit display that is easy to read even while wearing laser safety glasses.

High Energy Broadly Tunable OPA

FEATURES

- Ultra-wide spectral range from 193 to 16000 nm
- High peak power (>50 MW) ideal for non-linear spectroscopy applications
- Narrow linewidth <6 cm⁻¹ (for UV < 9 cm⁻¹)
- Motorized hands-free tuning in 193–2300 nm or 2300–16000 nm range
- ▶ PC control
- Remote control via keypad

APPLICATIONS

- Nonlinear spectroscopy: vibrational-SFG, surface-SH, Z-scan
- Pump-probe experiments
- ▶ Laser-induced fluorescence (LIF)
- Other laser spectroscopy applications

Available models

Model	Features
PG401	Model has a tuning range from 420 to 2300 nm and is optimized for providing highest pulse energy in the visible part of the spectrum. The wide tuning range makes PG401 units suitable for many spectroscopy application.
PG501-DFG	Model has a tuning range from 2300 to 16000 nm. The PG501-DFG1 model is the optimal choice for vibrational-SFG spectroscopy setups.

Picosecond Lasers

Other Ekspla Products



PGx01 SERIES

SPECIFICATIONS ¹⁾

Model	PG401	PG401-SH	PG401-DUV	PG501-DFG1	PG501-DFG2		
Tuning range							
DUV	-	-	193–209.95 nm	-	-		
SH	-	210–340, 370–419 nm		_			
Signal	420 – 680 nm		1	_			
Idler	740 – 2300 nm			-			
DFG				2300-10000 nm	2300-16000 nm		
Output pulse energy ²⁾	> 1000 µJ at 450 nm	> 100 µJ at 300 nm	> 50 µJ at 200 nm	> 250 µJ at 3700 nm, > 40 µJ at 10000 nm	> 250 µJ at 3700 nm, > 80 µJ at 10000 nm		
Linewidth	< 6 cm ⁻¹	< 9	cm⁻¹	< 6	cm ⁻¹		
Max pulse repetition rate			50 Hz				
Scanning step							
Signal	0.1 nm –						
Idler	1 nm			-			
Typical beam size ³⁾	~4 mm	~3	mm	~9	mm		
Beam divergence ⁴⁾		< 2 mrad		-	-		
Beam polarization	-	ver	tical	horiz	ontal		
Signal	horizontal			_			
Idler	horizontal						
Typical pulse duration	~20 ps	~15	5 ps	~20 ps			
PUMP LASER REQUIREMENTS							
Pump energy							
at 355 nm		10 mJ		-			
at 532 nm		-		10 mJ			
at 1064 nm	-	-	2 mJ	6 mJ	15 mJ		
Recommended pump source ⁵⁾		PL2231-50-TH, PL2251A-TH		PL2231-50-TH, PL2251A-TH	PL2231A-50-SH, PL2251B-SH		
Beam divergence			< 0.5 mra	d			
Beam profile		homogeneo	us, without hot spo	ots, Gaussian fit >90 %			
Pulse duration ⁶⁾			29 ± 4 ps				
PHYSICAL CHARACTERISTICS							
Size (W x L x H)	456 × 633 × 244 mm		456 × 10)31 × 249 ± 3 mm			
OPERATING REQUIREMENTS							
Room temperature			15 – 30 °C				
Power requirements		100 –	240 V AC single ph	ase, 47 – 63 Hz			
Power consumption			< 100 W				
 Due to continuous improvement, all specifications are subject to change w notice. Parameters marked typical are specifications. They are indications of t performance and will vary with each u manufacture. Unless stated otherwise, specifications are measured at 450 nm PG401 units, 3000 nm for PG501 units 300 nm for PG401SH units and for bas without options. See tuning curves for typical pulse energies available, please contact Ekspla for modetails. 	 Beam diameter is measured at the 1/e² level. * Full angle measured at the FWHM point. Full angle measured at the FWHM point. If a pump laser other than PL2250 or PL2230 is used, measured beam profile data should be presented when ordering. * Should be specified if non-EKSPLA pump laser is used. 				DANGER WELCHD OR KINDROUER TO DIRECT RELICTED OR SKITTERE RADUATION THE CONTRACT ON SKITTER RADUATION ON SKITTER RADUATION ON SKITTER RADUATION ON SK		

REV. 20220214 SAVANORIU AV. 237, LT-02300 VILNIUS, LITHUANIA TEL+370 5 2649629 E-MAIL SALES@EKSPLA.COM WWW.EKSPLA.COM



CUSTOMIZED FOR SPECIFIC REQUIREMENTS

Please note that these products are custom solutions tailored for specific applications or specific requirements.

Interested? Tell us more about your needs and we will be happy to provide you with tailored solution.

PG401-DFG1 provides:

51.0

- The broadest hands-free tuning range from 420 to 10000 nm
- It can be further extended up to 16000 nm with -DFG2 option. It should be noted, that for the 8000 – 16000 nm range a different nonlinear crystal is used, and exchange of the crystals needs to be done manually

Wavenumber, 10³ × cm⁻

49.0

48.0

50.0

PG402 features:

 Gap-free tuning range 410 – 709, 710 – 2300 nm

PGx01 SERIES

▶ Linewidth < 18 cm⁻¹

TUNING CURVES



Fig 1. Typical PG401 model tuning curve *Pump energy: 10 mJ at 355 nm*

Note: The energy tuning curves are affected by air absorption due narrow linewidth. These pictures present pulse energies where air absorption is negligible.

ON OPTICAL TABLE







Fig 4. Typical PG501-DFG1 tuning curve in 2300–10000 nm range *Pump energy: 7 mJ at 1064 nm*



Fig 3. Typical PG401-SH model tuning curve. Pump energy: 10 mJ at 355 nm



Fig 5. Typical PG501-DFG2 tuning curve in 8000–16000 nm range *Pump energy: 15 mJ at 1064 nm*



RECOMMENDED UNITS ARRANGEMENT

Fig 6. Arrangement of pump laser and PGx01 unit on optical table



Fig 7. Recommended arrangement of pump laser and PGx01-DFGx unit on optical table



PGx01 SERIES

OUTLINE DRAWINGS





OUTPUTS PORTS

Model	L, mm	a, mm	b, mm	c, mm	Port 1	Port 2
PG401	633	380	×	×	420–680 nm, 740–2300 nm	_
PG401-SH	838	380	×	×	210–340 nm, 370–419.9 nm, 420–680 nm, 740–2300 nm	_
PG401-SH/DUV	1026	380	250	50	210–340 nm, 370–419 nm, 420–680 nm, 740–2300 nm	192–209.95 nm

Fig 8. PG401 external dimensions



Fig 9. PG501 external dimensions

ORDERING INFORMATION

output port 2.

PG401-DUV

PG401-	DUV	PG501-DFG1
Model PG4xx → 355 nm pump	Optional tuning range extension DUV \rightarrow 193–209.95 nm SH \rightarrow 210-340 nm &	Model PG5xx → 532 nm pump
01 → travelling wave, narrowed linewidth 02 → travelling wave, net parrowed	370-420 nm Custom products, tailored for	01 → travelling wave, narrowed linewidth
not narrowed 11 → synchronous pumping, narrowed	specific applications. Inquire for other specifications.	Tuning range DFG1 → 2300–10000 nm;
	DFG1 → 2300–10000 nm; >250 µJ at 3700 nm DFG2 → 2300–16000 nm	>250 µJ at 3700 nm DFG2 → 2300−16000 nm

Note: Laser must be connected to the mains electricity all the time. If there will be no mains electricity for longer that 1 hour then laser (system) needs warm up for a few hours before switching on. Femtosecond Lasers



PICOSECOND TUNABLE SYSTEMS

PGx01 • PGx11 • PT277-XIR • PT277 • PT403

PGx11 SERIES



PGx11 series optical parametric devices employ advanced design concepts in order to produce broadly tunable picosecond pulses with nearly Fourier-transform limited linewidth and low divergence. High brightness output beam makes the PGx11 series units an excellent choice for advanced spectroscopy applications.

Optical layout of PGx11 units consists of Synchronously pumped Optical Parametric Oscillator (SOPO) and Optical Parametric Amplifier (OPA). SOPO is pumped by a train of pulses at approx. 87 MHz pulse repetition rate. The output from SOPO consists of a train of pulses

Features

spectroscopy.

series lasers with a 50 Hz pulse repetition rate.

with excellent spatial and spectral characteristics, determined by the SOPO cavity parameters.

OPA is pumped by a single pulse temporally overlapped with SOPO output. After amplification at SOPO resonating wavelength, the PGx11 output represents a high intensity single pulse on top of a low-intensity train, while in all other spectral ranges (idler for PG411 and PG711, signal for PG511, also DFG stages) only a single high intensity pulse is present.

Three models designed for pumping by up to the 3rd harmonic of Nd:YAG laser are available.

Model has a tuning range from 410 to 2300 nm and is optimized for providing highest pulse energy in the

visible part of the spectrum. When combined with an optional Second Harmonic Generator (SHG) and Sum

Model has a tuning range 2300–10000 nm. PG411 and PG511 models are designed to be pumped by PL2230

Model has 1 kHz pulse repetition rate and uses DPSS mode-locked laser of the PL2210 series for pumping. When pumped with pulses of 90 ps duration, linewidths of less than 1 cm⁻¹ were measured in the spectral

range up to 16 µm, which makes this device an excellent choice for time-resolved or nonlinear infrared

Frequency Generator (-DUV), it offers the widest possible tuning range - from 193 to 2300 nm.

Transform Limited Broadly Tunable Picosecond OPA

FEATURES

- > 2 cm⁻¹ or **1 cm⁻¹** linewidth
- High brightness picosecond pulses at 50 Hz or at up to 1 kHz pulse repetition rate
- Nearly Fourier-transform limited linewidth
- Low divergence <2 mrad</p>
- ▶ Hands-free wavelength tuning
- Tuning range from 193 nm to 16000 nm
- PC control
- Remote control via keypad

APPLICATIONS

- Time resolved pump-probe spectroscopy
- Laser-induced fluorescence
- Infrared spectroscopy
- Nonlinear spectroscopy: vibrational-SFG, surface-SH, Z-scan, pump probe
- Other laser spectroscopy applications

32



Available models

Model

PG411

PG511

PG711

Femtosecond Lasers



Microprocessor based control system provides automatic positioning of relevant components, allowing hands free operation. Nonlinear crystals, diffraction grating and filters are rotated by ultra-precise stepper motors in microstepping mode, with excellent reproducibility. Precise nonlinear crystal temperature stabilization ensures long-term stability of generated wavelength and output power.

For customer convenience the system can be controlled through USB (VCP, ASCII commands), RS232 (ASCII commands), LAN (REST API) or RS232 (ASCII commands), LAN (REST API) depending on the system configuration or a remote control pad. Both options allow easy control of system settings.

Available standard models are summarized in a table below. Please inquire for custom-built versions.

SPECIFICATIONS ¹⁾

Model	PG411 PG411-SH PG411-SH-D		PG411-SH-DUV	PG511-DFG	PG711	PG711-DFG
Output wavelength tuni	ng range					
SH, DUV	-	210–410 nm	193–410 nm		-	
Signal		410–709 nm	้า	-	1550-	2020 nm
Idler		710–2300 nr	n	-	2250-	-3350 nm
DFG	_			2300–10000 nm	-	3350–16000 nm
DFG2 (up to 16000 nm)	_			inquire	-	
Output pulse energy ²⁾						
SH, DUV	-	100 µJ 3)	50 µJ 3)	-		
Signal		700 µJ		-	5	00 μJ
Idler 4)		250 µJ		-	1(Lμ 0C
DFG	_			> 200 µJ at 3700 nm, > 40 µJ at 10000 nm	_	20 µJ 5)
Max pulse repetition rate	50 Hz			50 Hz	1000 Hz	
Linewidth	< 3 cm ^{-1 6)}			< 2 cm ⁻¹	<	1 cm ⁻¹
Linewidth Idler		< 5 cm ^{-1 6)}			-	
Typical pulse duration ⁷⁾		~20 ps		~20 ps	~	70 ps
Scanning step						
SH, DUV	-	0	.01 nm		-	
Signal				0.1 nm		
ldler				1 nm		
DFG			-			1 nm
Typical beam diameter ⁸⁾		~ 4 mm		~ 9 mm	~	3 mm
Beam divergence 9)				< 2 mrad		
Beam polarization ⁹⁾						
SH, DUV	-	v	rertical		-	
Signal		horizontal		vertical horizontal		
ldler		vertical		horizontal	Ve	ertical
DFG	_			horizontal	-	horizontal

¹⁾ Due to continuous improvement, all specifications are subject to change without notice. Parameters marked typical are not specifications. They are indications of typical performance and will vary with each unit we manufacture. Unless stated otherwise, all specifications are measured at 450 nm for PG411 units, 800 nm for PG511 units, and 1620 nm for PG711 units and for basic system without options.

- ²⁾ Pulse energies are specified at selected wavelengths. See typical tuning curves for pulse energies at other wavelengths.
- ³⁾ Measured at 280 nm for SH and 200 nm for DUV.
- ⁴⁾ Measured at 1000 nm for PG411 units, 1620 nm for PG511, and 3000 nm for PG711 units.

⁵⁾ Measured at 10000 nm.

- $^{6)}$ Linewidth for signal (409 710 nm) < 3 cm $^{1},$ linewidth for idler and SH-DUV (710 2300 nm and 193 409 nm) < 5 cm $^{1}.$
- ⁷⁾ Estimated FWHM assuming pump pulse duration 30 ps at 1064 nm for PG411 and PG511 units, and 90 ps at 1064 nm for PG711 units.
- 8 Beam diameter is measured at 1/e² level and can vary depending on the pump pulse energy.
- ⁹⁾ Full angle measured at FWHM level.





PGx11 SERIES

SPECIFICATIONS ¹⁾

Model	PG411	PG411-SH	PG411-SH-DUV	PG511-DFG	PG711	PG711-DFG		
PUMP LASER REQUIR	REMENTS							
Recommended pump source	PL223	PL2231 + APL2100-TRAIN-H411 PL2231 + H500- APL2100-TRAIN PL2211A TR						
Min. pump energy or po	ower ¹⁰⁾							
at 1064 nm	-		2 mJ	(10 mJ)	Eml	at 1 kHz		
at 532 nm		-		5 mJ (8 mJ)	נוח כ			
at 355 nm		5 mJ (10 m.	J)		-			
Pulse duration ¹¹⁾			30 ps		ç	00 ps		
Bream polarization at pump wavelength		vertical horizontal						
Beam size 12)		7 mm 2.5 mm						
Beam divergence		< 0.5 mrad						
Beam profile			homogeneo	us, without hot spots				
PHYSICAL CHARACTE	RISTICS							
Size (W \times L \times H)	456 × 1026 × 244 mm	456 × 1026 × 244 mm 456 × 1226 × 244 mm PL2231: 456 × 1026 × 244 mm 456 × 244 mm 456 × 1026 × 244 mm 456 × 1026 × 244 mm 456 × 24						
OPERATING REQUIRE	MENTS							
Room temperature				15–30 °C				
Room temperature stability		± 2 °C						
Power requirements	100–240 V single phase, 47–63 Hz							
Power consumption		< 300 W						

¹⁰⁾ The first number represents pulse train energy or power, while the value in brackets represents single pulse energy. ¹¹⁾ At FWHM level. Inquire for other available pulse duration options.

 $^{\scriptscriptstyle 12)}\,$ Beam diameter measured at 1/e² level.

RECOMMENDED UNITS ARRANGEMENT ON OPTICAL TABLE



Fig 1. Arrangement of pump laser and PGx11 unit on optical table

Picosecond Lasers

34

PGx11 SERIES

TUNING CURVES













Fig 5. Typical PG711-DFG model tuning curve. *Pump energy: 2.5 mJ at 1064 nm, 1 kHz repetition rate*



Note: The energy tuning curves are affected by air absorption due narrow linewidth. These pictures present pulse energies where air absorption is negligible.





PGx11 SERIES

OUTLINE DRAWINGS





Output panel

1

Т

Fig 7. PG411 model external dimensions



Model	L, mm	a, mm	b, mm	c, mm	Port 1	Port 2	Port 3
PG411	1026	×	411	×	420–709 nm, 710–2300 nm	420–709 nm, 710–2300 nm	_
PG411-SH	1226	×	411	×	420–709 nm, 710–2300 nm	210–419 nm, 420–709 nm, 710–2300 nm	-
PG411-SH/DUV	1226	235	411	331	420–709 nm, 710–2300 nm	210–419 nm, 420–709 nm, 710–2300 nm	192–209.95 nm

ORDERING INFORMATION

Note: Laser must be					
connected to the mains					
electricity all the time. If there					
will be no mains electricity for					
longer that 1 hour then laser					
(system) needs warm up for a					
few hours before switching on.					

PGx11-SH

Model	Optional tuning range extension
PG411 \rightarrow ps 355 nm pump	SH (PG411) → 210-420 nm
PG511 → ps 532 nm pump	SH/DUV (PG411) → 193-420 nm
PG711 → ps 1064 nm pump	DFG (PG511) → 2300–10000 nm
	DFG (PG711) \rightarrow 3350–16000 nm



PICOSECOND TUNABLE SYSTEMS

PGx01 • PGx11 • PT277-XIR • PT277 • PT403

PT277-XIR SERIES



PT277 series laser systems integrate a picosecond optical parametric oscillator and DPSS pump laser into a single compact housing. Mounting the components into one frame provides a cost-effective and robust solution. It makes laser installation shorter and improves long-term stability, reduces maintenance costs. The tuning range for the model PT277-XIR 1400 – 2000, 2280 – 4450, 4800 – 16000 nm (7000 – 5000, 4385 – 2245, 2080 – 625 cm⁻¹) with linewidth <5 cm⁻¹ in the full tuning range. The fast wavelength tuning is based on the microprocessor-control and wavelength tuning is fully automatic. The wavelength tuning elements are mounted on precise closed loop micro stepping motors. The temperatures of the non linear crystals is controlled by a precise thermocontrollers. For customer convenience the system can be controlled via keypad and/or any controller running on any OS using REST API commands. Variety of interfaces USB, RS232, LAN, WLAN (optionally) ensures easy control and integration with other equipment.

Single Housing MIR (mid infrared range) Tunable Picosecond Laser

FEATURES

- 1400 16000 nm (7000 – 625 cm⁻¹) tuning range
- Linewidth <5 cm⁻¹ in the full tuning range
- Nearly diffraction limited divergence
- Remote control via keypad and/or any controller running on any OS using REST API commands
- PC control via USB (virtual com port), RS232, LAN

APPLICATIONS

▶ Infrared spectroscopy

TUNING CURVES

Fig 1. Typical output PT277-XIR tuning curve. *The actual tuning curve might differ from presented here.*





Picosecond Lasers

Picosecond Tunable Systems

SPECIFICATIONS

Femtosecond Lasers

Picosecond Lasers

Picosecond Tunable Systems

Nanosecond Lasers

Model	PT277 - XIR		
SYNCHRONOUSLY PUMPED			
Tuning range			
Signal	1400 – 2000 nm (7000 – 5000 cm ⁻¹)		
Idler	2280 – 4450 nm (4385 – 2245 cm ⁻¹)		
XIR	4800 – 16000 nm (2080 – 625 cm ⁻¹)		
Linewidth	< 5 cm ⁻¹		
Output power 2)			
OPO @ 4000 nm (2500 cm ⁻¹)	> 150 mW		
OPO @ 12500 nm (800 cm ⁻¹)	> 10 mW		
Pulse repetition rate	87 MHz		
Pulse duration (pump laser)	~8 ps		
Scanning step			
Signal	0.1 nm		
ldler, XIR	1 nm		
Polarization	vertical		
Typical beam diameter ^{3) 4)}	~3 mm		
Typical beam divergence 5)	< 2 mrad (for signal)		
Beam pointing stability	< 50 µrad rms @ 1596 nm		
AOM modulation	70 kHz – 2 MHz		
Fast spectral scan, scan speed (for spectral range)			
Idler	< 2 s		
XIR	< 1 s		
PHYSICAL CHARACTERISTICS	5		
Unit size (W×L×H)	320 × 766 × 241 mm		
Power supply size (W×L×H)	483 × 140 × 390 mm		
Umbilical length	2.5 m		
OPERATING REQUIREMENTS			
Cooling	water-air (by provided chiller)		
Room temperature	22 ± 2 °C		
Room temperature stability	±1°C		
Relative humidity	20 – 80 % (noncondensing)		
Power requirements	100 - 240 VAC (-10% / +5%), single phase, 50/60 Hz		
Power consumption	< 1 kVA		
Cleanness of the room	not worse than ISO Class 9		



- 1) All specifications are subject to change without notice. The parameters marked typical are not specifications. They are indications of typical performance and will vary with each unit we manufacture. Unless stated otherwise all specifications are measured at 1596 nm.
- ²⁾ Output powers are specified at selected wavelengths. See typical tuning curves for power at other wavelengths.
- ³⁾ Measured at 3000 nm.
- Beam diameter at the 1/e² level and can vary 4) depending on the pump pulse energy.
- ⁵⁾ Full angle measured at the FWHM level.

Note: Laser must be connected to the mains electricity all the time. If there will be no mains electricity for longer that 1 hour then laser (system) needs warm up for a few hours before switching on. The laser and auxiliary units must be settled in such a place void of dust and aerosols. It is advisable to operate the laser in air conditioned room, provided that the laser is placed at a distance from air conditioning outlets. The laser should be positioned on a solid worktable. Access from both sides should be ensured. Intensive sources of vibration should be avoided, like railways station etc nearby laboratory.







OUTLINE DRAWINGS

PICOSECOND TUNABLE SYSTEMS

PGx01 • PGx11 • PT277-XIR • PT277 • PT403

PT277 SERIES



PT277 series laser systems integrate a picosecond optical parametric oscillator and DPSS pump laser into a single compact housing. Mounting the components into one frame provides a cost-effective and robust solution with improved long-term stability and reduced maintenance costs.

The tuning range is for the model PT277 1400 – 2050 and 2200 to 4450 nm with nearly Fourier transform limited linewidth.

The microprocessor-controlled wavelength tuning is fully automatic. The wavelength controlling elements are mounted on precise micro-stepping motors. The temperature of the non-linear crystal is controlled by a precise thermocontroller with a bidirectional Peltier element, resulting in the fast tuning of crystal temperature. For customer convenience the system can be controlled through USB (VCP, ASCII commands), RS232 (ASCII commands), LAN (REST API) or RS232 (ASCII commands), LAN (REST API) depending on the system configuration or a remote control pad. Both options allow easy control of system settings.

Single Housing NIR-IR Range Tunable Picosecond Laser

FEATURES

- ▶ 1400–4450 nm tuning range
- Nearly Fourier transform-limited linewidth
- Nearly diffraction limited divergence
- Output wavelength monitoring (optional)
- ▶ PC control

APPLICATIONS

- ▶ Infrared microscopy
- ▶ Infrared spectroscopy
- Near field spectroscopy

TUNING CURVES

Fig 1. Typical output power of PT277 tunable laser. The power is shown only at the wavelengths where ambient air absorption is negliglible





PT277 SERIES

SPECIFICATIONS ¹⁾

Model	PT277		
Pulse repetition rate ²⁾	87 MHz		
Tuning range			
Signal	1400 – 2050 nm		
ldler	2200 – 4450 nm		
Output power ³⁾			
OPO ⁴⁾	> 500 mW		
Linewidth 4)	< 2.5 cm ⁻¹		
Typical pulse duration ^{4) 5)}	70 ps		
Scanning step			
Signal	0.1 nm		
ldler	0.1 nm		
Polarization			
Signal beam	horizontal		
Idler beam	horizontal		
Typical beam diameter ^{4) 6)}	~2 mm		
Typical beam diameter, Idler 4) 6)	~5 mm		
Typical beam divergence ^{4) 7)}	< 2 mrad		
PHYSICAL CHARACTERISTICS			
Unit size (W \times L \times H)	370 × 800 × 260 mm		
Power supply size (W \times L \times H)	520 × 500 × 290 mm		
Umbilical length	2 m		
OPERATING REQUIREMENTS			
Cooling	water-air		
Room temperature	22 ± 2 °C		
Relative humidity	20 – 80 % (noncondensing)		

Power consumption

Power requirements

- Due to continuous improvement, all specifications are subject to change without notice. Parameters marked 'typical' are indications of typical performance (not specifications) and will vary with each unit we manufacture. Unless stated otherwise, all specifications are measured at 1064 nm and for basic system without options.
- ²⁾ Inquire for custom pulse repetition rates.
- ³⁾ Output powers are specified at selected wavelengths. See typical tuning curves for power at other wavelengths.

⁴⁾ Measured at 1600 nm for PT277 model at signal range.

- ⁵⁾ Pulse duration can vary depending on wavelength and pump energy.
- ⁶⁾ Beam diameter at the 1/e² level and can vary depending on the pump pulse energy.
- 7) Full angle measured at the FWHM level.



OUTLINE DRAWINGS

Note: Laser must be connected to the mains electricity all the time. If there will be no mains electricity for longer that 1 hour then laser (system) needs warm up for a few hours before switching on.



Fig 2. External dimensions of PT277 series lasers



100 – 240 V AC, single phase 50/60 Hz

< 1 kVA



PGx01 • PGx11 • PT277-XIR • PT277 • PT403

PT403 SERIES



PT403 series laser systems integrate a picosecond 1 kHz repetition rate DPSS pump laser and optical parametric generator into a single housing. New picosecond tunable wavelength laser system provide from 210 to 2300 nm from the one box.

Unlike other solutions in the market, offering laser and OPO in different units, new approach features pump laser and OPO integrated into one unit. That delivers almost twice smaller footprint, shorter installation, better stability and other substantial benefits for user.

All-in one-box solution features all components placed into one compact housing. It means better overall stability because all potential causes for misalignment between separate units of pump laser and optical parametric generator are eliminated.

To ensure reliability industry and market tested solutions were employed during the build-up of PT403.

Pump laser is based on industry "gold standard" diode pumped Ekspla PL2210 series picosecond mode-locked laser. Improved output parameters and reduced maintenance costs are achieved by employing diode-pumped-only technology. Optical parametric generator is based on PGx03 picosecond optical parametric amplifier systems. Fully automatized and microprocessor based control system ensures hands free precise wavelength tuning.

PT403 was built without sacrificing any parameters or reliability. The optical design is optimized to produce low divergence beams with moderate linewidth (typically < 9 cm⁻¹) at approximately 20 ps pulse duration. Featuring 1 kHz repetition rate PT403 tuneable laser is versatile cost-efficient tool for scientists researching various kind of disciplines like time resolved fluorescence, pump-probe spectroscopy, laserinduced fluorescence, Infrared spectroscopy and other aplications.

Simple and convenient laser control

For customer convenience the laser can be operated from master device or personal computer through USB (VCP, ASCII commands), RS232 (ASCII commands), LAN (REST API) or RS232 (ASCII commands), LAN (REST API) depending on the system configuration or from remote control pad with backlit display that is easy to read even while wearing laser safety glasses.

Tunable Wavelength Picosecond Laser

FEATURES

- Tuning range: 210 2300 nm
- Motorized hands-free tuning
- High pulse energy at 1 kHz rates
- Diode pumped solid state design
- Narrow linewidth < 9 cm⁻¹
- Remote control via keypad
- ▶ PC control
- Optional streak camera triggering pulse with < 10 ps rms jitter
- ▶ Turn-key operation
- Air cooled external water supply is not required
- Low maintenance costs

APPLICATIONS

- Time resolved fluorescence (including streak camera measurements), pump-probe spectroscopy
- Laser-induced fluorescence
- ► Infrared spectroscopy
- Nonlinear spectroscopy: surface-SH, Z-scan
- Other spectroscopic and nonlinear optics applications

BENEFITS

- Better long term stability (compared with layout where laser and OPO are in different units)
- Higher safety all beams are in the box
- Shorter installation time
- Almost twice smaller footprint



Picosecond Tunable Systems

Nanosecond Lasers

Nanosecond Tunable Lasers

PT403 SERIES

SPECIFICATIONS ¹⁾

Model	PT403	PT403-SH				
OPA SPECIFICATIONS						
Output wavelength tuning range						
SH	-	210 – 409 nm				
Signal	410 – 709 nm					
Idler	710 – 2300 nm					
Output pulse energy ²⁾						
SH ³⁾	_	15 µJ				
Signal 4)	لµ 08	· · · · ·				
Idler ⁵⁾	25 µJ					
Pulse repetition rate	1000 Hz					
Linewidth	< 9 cm ⁻¹	< 12 cm ⁻¹				
Typical pulse duration ⁶⁾	~ 20 ps					
Scanning step						
SH	_	0.05 nm				
Signal	0.1 nm					
Idler	1 nm					
Typical beam size 7	~ 2 mm					
Beam divergence ⁸⁾	< 2 mrad					
Beam pointing stability	≤ 100 µrad r	< 100 µrad rms				
Beam polarization	•					
SH	_	horizontal				
Signal	horizontal					
Idler	vertical	vertical				
Optical pulse jitter						
Internal triggering regime ⁹⁾	< 50 ps (StDev) in respect to	o TRIG1 OUT pulse				
External triggering regime	~ 3 ns (StDev) in respect t	o SYNC IN pulse				
TRIG1 OUT pulse delay ¹⁰⁾	-400 150 ns					
ODEDATING REQUIREMENTS						
Room temperature	22 + 2 °C					
Relative humidity	20 – 80% (non-cor	ndensina)				
Power requirements	100 – 240 V single phase 47 – 63 Hz					
Power consumption	<pre></pre>					
Water service	air cooled					
Cleanness of the room	not worse than ISO Class 9					
¹⁰ Due to continuous improvement, all specifications are subject to change without notice. Parameters marked typical are not specifications. They are indications of typical performance and will vary with each unit we manufacture. Unless stated otherwise, all specifications are measured at 450 nm for PT403 units for basic system without options.	 ³⁾ Measured at 260 nm. ⁴⁾ Measured at 450 nm. ⁵⁾ Measured at 1000 nm. ⁶⁾ Estimated assuming 30 ps at 1064 nm pump pulse. Puls duration varies depending on wavelength and pump er ⁷⁾ Beam diameter at the 1/e² level. Can vary depending o wavelength 	See wished and/org invisible LASER RADIATION with the first of the contract of				

²⁾ Pulse energies are specified at selected wavelengths. See typical tuning curves for pulse energies at other wavelengths.

⁸⁾ Beam divergence measured at FWHM.

 $_{\rm 9)}\,$ < 10 ps jitter is provided with PRETRIG option.

Communication module interfaces

Interface	Description				
USB	virtual serial port, ASCII commands				
RS232	ASCII commands				
LAN	REST API				
WLAN	REST API				



Picosecond Lasers

Other Ekspla Products

42

PT403 SERIES

DESIGN

The units can be divided into several functional parts:

- 1. 1 kHz repetition rate DPSS pump laser,
- 2. Optical parametric generator (OPG),
- 3. Electronic control unit.



Fig 1. PT403 unit

PT403 series laser systems integrate a picosecond 1 kHz repetition rate DPSS pump laser and optical parametric generator into a single housing. As pump laser is used PL2210 series diode-pumped, air-cooled, mode-locked Nd:YAG laser. Picosecond tunable wavelength laser system provide from 210 to 2300 nm from the single optical unit.

OPTIONS

Option SF

Energy increasing in 300 – 409 nm range by sum-frequency generation.

> 20 μJ @ 340 nm. Pulse energies are \sim 10 % lower in comparison to the system without SE option. See table below for pulse energy specifications:

	without Sr option. See ta	Die Delow for pulse	energy s	pecifications.
	Model ¹⁾	PT403		PT403-SH
	SH ²⁾	-		> 13 µJ
	Signal ³⁾		> 7	lμ 0
	Idler ⁴⁾		> 2	2 µJ
	 Due to continuous improvem subject to change without no typical are not specifications. typical performance and will w manufacture. Options -H, -2H, -3H 	ent, all specifications are tice. Parameters marked They are indications of rary with each unit we	 e ² Measured at 260 nm. d ³ Measured at 450 nm. 4 Measured at 1000 nm. 	
	1064 nm or 532 nm, or 35 – H output energy 0.7 mJ – 2H output energy 0.3 m – 3H output energy 0.3 m	55 nm outputs ^{າງ 2)} ; ນ; ນ.	¹⁾ Ou ²⁾ Inq with	tputs are not simultaneous. uire for outputs simultaneously n PG.
CUSTOMIZED FOR SPECIFIC REQUIREMENTS	Please note that these pr custom solutions tailored applications or specific re	oducts are I for specific equirements.	Intereste your need provide y	d? Tell us more about Is and we will be happy to ou with tailored solution.

PT503 FEATURES

▶ The higher pulse energy in the near-IR spectral range ▶ Tuning range from 700 to 2200 nm

Femtosecond Lasers



PT403 SERIES

TUNING CURVES



Fig 2. Typical PT403 tuning curves in signal (420 – 709 nm) and idler (710 – 2300 nm) ranges.

Note: The energy tuning curves are affected by air absorption due narrow linewidth. These pictures present pulse energies where air absorption is negligible.

OUTLINE DRAWINGS



Fig 3. PT403 series laser head typical outline drawing



Fig 4. Compared with layout where laser and OPO are in different units, PT403 features almost twice smaller footprint

www.gmp.ch

s 17	CH-1020 Renens	Tél. 021 633 21 21	Fax. 021 633 21 29	633 21 29 info@gmp.ch			
a	CH-8117 Fällanden	Tel. 044 825 34 00	Fax. 044 825 34 01 info@gmp.ch		Port 2	Port 3	Port 4
	PT403	1241	1064	1 / 532 nm	-	355 nm	410 – 2300 nm
	PT403-SH/SF	1441	1064	1 / 532 nm	210 – 2300 nm	355 nm	410 – 2300 nm

Note: Laser must be connected to the mains electricity all the time. If there will be no mains electricity for longer that 1 hour then laser (system) needs warm up for a few hours before switching on.

Picosecond Lasers

e: Avenue des Baumette ch: Dübendorfstrasse 11





www.gmp.ch

GMP SA GMP SA

Main office: Avenue des Baumettes 17 Büro Zürich: Dübendorfstrasse 11a

CH-1020 Renens CH-8117 Fällanden

Tél. 021 633 21 21 Tel. 044 825 34 00

Fax. 021 633 21 29 Fax. 044 825 34 01 info@gmp.ch info@gmp.ch