Nitrogen Laser

NL100 — 337 nm nitrogen laser



- 337 nm wavelength
- 170 µJ pulse energy
- Internal or external triggering to 20 Hz
- Replaceable laser cartridge
- No mirror alignment necessary
- TTL level sync output pulse (opt.)
- Fully compatible with VSL-337i OEM

NL100 Nitrogen Laser

The NL100 Nitrogen Laser is ideal for fluorescence measurements, MALDI-TOF mass spectrometers, and other pulsed UV radiation experiments. It provides 3.5 ns pulses at 337 nm (UV), with repetition rates up to 20 Hz. The pulse energy is $170\,\mu J,$ which results in a peak power of $45\,kW$ and an average power of 3 mW.

The NL100 can be triggered internally or externally at rates up to 20 Hz. It can also provide a sync output pulse (optional) derived from the laser pulse for experiments or systems where sub-nanosecond accuracy is critical. The user also has the option of running the laser system in command charge mode.

The NL100 uses a replaceable, sealed laser cartridge which includes the high voltage storage capacitors, switching element, and laser tube. The cartridge is warranted to maintain at least 70% of its energy for twenty million pulses or two years, whichever occurs first.

No mirror alignment is ever necessary in the NL100, as the laser optics are mounted on the plasma tube and aligned at the factory. The NL100 also includes all safety features necessary to comply with the U.S. laser safety standards contained in 21 CFR 1040.10.

The NL100 is fully compatible with the Spectra-Physics model VSL-337i OEM Nitrogen Laser (part #337999-04). The mechanical design of the NL100 includes all of the mounting and alignment features of the VSL-337i, making the NL100 a straightforward retrofit in many pre-existing systems.

www.gmp.ch

GMP SA Main office: Avenue des Baumettes 17 GMP SA Büro Zürich: Dübendorfstrasse 11a

CH -1020 Renens CH-8117 Fällanden

Tel. 044 825 34 00

Tél. 021 633 21 21

info@gmp.ch Fax. 044 825 34 01 info@gmp.ch

Fax. 021 633 21 29

Beam Characteristics

Wavelength Spectral bandwidth Pulse width (FWHM) Pulse energy Energy stability (pulse to pulse) Peak power Average power Beam size Beam divergence (full angle)

Triggering

Repetition rate

External trigger input Internal trigger generator Sync output pulse

General

Power requirements

Power consumption Key switch Interlock switch Dimensions

Weight Warranty

lion

337.1 nm 0.1 nm <3.5 ns 170 μJ 3% std. dev. (at 10 Hz) 45 kW 3 mW (at 20 Hz) 3×7 mm 5×8 mrad

0 to 20 Hz (external trigger) 1 to 20 Hz (internal trigger) TTL (opto-isolated) 1 to 20 Hz (adjustable) TTL level (opt.)

+24 VDC, 1.5 A (average) at 20 Hz., 3 A (peak) 36 W (20 Hz operation) On/off Built-in 3.75"×3.75"×11" (WHD) (9.5 cm×9.5 cm×27.9 cm) 7.5 lbs., 3.4 kg Cartridge is warranted to maintain at least 70% of its energy for twenty milpulses or two years, whichever occurs first.



NL100 front panel



NL100 rear panel

Ordering Information

NL100	Nitrogen laser
O100SP	Optional sync-out pulse
O100RC	Replacement laser cartridge

www.gmp.ch

GMP SA Main office: Avenue des Baumettes 17 GMP SA Büro Zürich: Dübendorfstrasse 11a CH -1020 Renens CH-8117 Fällanden

 Renens
 Tél. 021 633 21 21

 Fällanden
 Tel. 044 825 34 00

Fax. 021 633 21 29 info@gmp.ch Fax. 044 825 34 01 info@gmp.ch