



Import, delivery and installation of optical tables



TMC all-steel optical tables and vibration isolation lab tables can weight several hundred kilos and must be handled and assembled with care and precision.




GMP offers a full installation service on demand.





From Factory to Laboratory a full service package from GMP



Package 1 Sales and delivery

- implementation study on site
- shipping, customs clearance
- import duties and insurance
- transport cost to the final destination loading dock
- warranty in Switzerland

calculated with the price of the table

Package 2 Optional full installation

- visit with transport specialist
- rental of crane if necessary
- floor protection
- transfer into the laboratory
- unpacking
- assembling by a technical expert
- connection to a compressed air supply
- control and adjusting horizontally

quote on request depending on location, size and weight

References

«L'installation de 4 tables optiques a été exécutée avec calme, efficacement et avec un grand professionnalisme par M. P.-A. Thonney, le technicien de la Maison GMP. Le personnel de GMP s'est montré très flexible et avenant lors des changements de dates dus à des retards dans la construction des Laboratoires. Le personnel, vendeurs inclus, est très compétent et très serviable»

J.-C.T., ETH Zentrum

«Vielen Dank für die problemlose Installation!»

M. H., ETH Hönggerberg

«Four anti-vibration tables were installed in our laboratory with professionalism and efficiency, as was the (very rare) technical tuning to maintain them at top performance»

Prof. R.B., Université de Lausanne